

# Boule Cal

**CAL****IVD**

Calibrator

| Lot      |  |
|----------|-----------------------------------------------------------------------------------|
| 22404-34 | 2024-07-14                                                                        |

Open vial stability 5 days <sup>1</sup>**Medonic**

| Parameter <sup>2</sup> | RBC 10 <sup>12</sup> /L | MCV fL | Plt 10 <sup>9</sup> /L | MPV fL | WBC 10 <sup>9</sup> /L | Hgb g/dL |
|------------------------|-------------------------|--------|------------------------|--------|------------------------|----------|
| Medonic M-series / M32 | 4.12                    | 74.2   | 238                    | 11.9   | 8.5                    | 11.0     |
| Range ± <sup>3</sup>   | 0.10                    | 2.0    | 15                     | 1.5    | 0.3                    | 0.2      |

**Swelab**

| Parameter <sup>2</sup>  | RBC 10 <sup>12</sup> /L | MCV fL | Plt 10 <sup>9</sup> /L | MPV fL | WBC 10 <sup>9</sup> /L | Hgb g/dL |
|-------------------------|-------------------------|--------|------------------------|--------|------------------------|----------|
| Swelab Alfa / Alfa Plus | 4.12                    | 75.0   | 243                    | 12.9   | 8.5                    | 11.0     |
| Range ± <sup>3</sup>    | 0.10                    | 2.0    | 15                     | 1.5    | 0.3                    | 0.2      |

|    | 1                                            |
|----|----------------------------------------------|
| bg | Трайност на отворен флакон                   |
| hr | Stabilnost otvorene bočice                   |
| cs | Stabilita otevřené ampulky                   |
| da | Åben rør stabilitet                          |
| et | Stabiilsus avatud viaalis                    |
| fr | Stabilité en flacon ouvert                   |
| de | Stabilität geöffneter Flaschen               |
| el | Διατήρηση δείγματος μετά την αποσφράγιση     |
| it | Stabilità della fiala aperta                 |
| lv | Atvērt flakonu stabilitāti                   |
| lt | Stabilumas atidarius buteliuką               |
| no | Åpen rør/glass stabilitet                    |
| pl | Trwałość po otwarciu fioleki                 |
| pt | Estabilidade após abertura do frasco         |
| ro | Termenul de valabilitate al fiolei desfăcută |
| ru | Стабильность вскрытого флакона               |
| sr | Proizvod stabilan nakon otvaranja pakovanja  |
| sk | Stabilita otvorenej injekčnej liekovky       |
| sl | Stabilnost odprte stekleničke                |
| es | Estabilidad de la cápsula abierta            |
| sv | Hållbarhet för öppnad flaska                 |
| tr | Açık şişe stabilitesi                        |

|               | 2          | 3                   |
|---------------|------------|---------------------|
| 5 дни         | Параметър  | Обхват              |
| 5 dni         | Parametar  | Raspon              |
| 5 dny         | Parametr   | Rozsah              |
| 5 dage        | Parameter  | Område              |
| 5 päevade arv | Parameeter | Vahemik             |
| 5 jour        | Paramètre  | Intervalle          |
| 5 tage        | Parameter  | Bereich             |
| 5 ημέρα       | Παράμετρος | Αναμενόμενο εύρος α |
| 5 giorni      | Parametro  | Intervallo          |
| 5 dienas      | Parametrs  | Diapazons           |
| 5 d.          | Parametras | Intervalas          |
| 5 dager       | Parameter  | Område              |
| 5 Liczba dni  | Kryterium  | Zakres              |
| 5 dias        | Parâmetro  | Intervalo           |
| 5 zile        | Parametru  | Interval            |
| 5 дней        | Параметр   | диапазон            |
| 5 dana        | Parametar  | Opseg               |
| 5 dní         | Parameter  | Rozsah              |
| 5 dni         | Parameter  | Območje             |
| 5 dias        | Parâmetro  | Intervalo           |
| 5 dagar       | Parameter  | Intervall           |
| 5 gün         | Parametre  | Aralığı             |

[www.boule.com](http://www.boule.com)

 **Boule**  
Boule Medical AB  
Domnarvsgatan 4  
SE-163 53 Spånga, Sweden

25123-6

1504055

## For Medonic M-series / M32 with barcode reader

### M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

### M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series

Lot: 22404-34

Prod.: Cal

ID: 2240434+

Exp.: 14-Jul-2024



## For Swelab Alfa / Alfa Plus with barcode reader

### Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

### Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa

Lot: 22404-34

Prod.: Cal

ID: 2240434+

Exp.: 14-Jul-2024

