

Certificate of Analysis

| | |
|---------------------------|--|
| Product | SWELAB ALFA COMBOPACK (1 L Lyse / 5 L Diluent) |
| Product No | 1504127 |
| Lot No | 2405-196 |
| Manufacturing Date | 2024-05-21 |
| Expiration Date | 2027-05-20 |

| Test | Specification Diluent | Specification Lyse |
|------------------------------|------------------------------|---------------------------|
| pH | 6,70 – 6,90 | 7,80 – 8,30 |
| Osmolality (mOsm/kg) | 294 – 312 | 295 – 330 |
| Conductivity (mS/cm) | 16,98 – 17,93 | 17,75 – 18,25 |
| CBC Performance * | Passes test | Passes test |
| Particles (PLT/WBC channels) | ≤ 5 | ≤ 0,2 |

* Performance test using calibrated Blood cell counter at Boule Medical reference laboratory. All parameters except mid cells are tested. Accepted reference 95 % to 105 % of expected control blood sample values.

Actual values meet specifications.

Date: 2024-05-28



Caroline Sparf
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

