

LOT: 12405-52
EXP: 2024-10-20
Open Vial Stability: 16 days

EXIGO		Mean	± Range	Expected Range	
WBC	10 ⁹ /L	8.6	0.6	8.0	9.2
Hct	%	28.9	2.5	26.4	31.4
MCV	fL	67.8	5.0	62.8	72.8
RBC	10 ¹² /L	4.26	0.18	4.08	4.44
Hgb	g/dL	10.2	0.4	9.8	10.6
MCH	pg	23.9	1.5	22.4	25.4
MCHC	g/dL	35.3	2.2	33.1	37.5
PLT	10 ⁹ /L	207	30	177	237

This assay sheet is also available on :

www.boule.com

For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12405-52
Prod.: Vet-Con Norm
ID: 1240552+
Exp.: 20-Oct-2024

