

Boule Con-5Diff A1

CONTROL

IVD

LOT



Tri-Level	72405-1K	2024/09/05
Low	72405-11	2024/09/05
Normal	72405-12	2024/09/05
High	72405-13	2024/09/05

Open vial stability: 14 days ¹

Medonic M51 / Swelab Lumi

Parameter ²	Unit ³	Low		Normal		High	
		Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵
WBC	10 ⁹ /L	3.40	0.50	8.08	1.00	18.75	2.50
Neu%	%	49.3	16.0	53.9	8.0	62.8	8.0
Lym%	%	40.0	21.0	31.3	8.0	20.0	12.0
Mon%	%	4.1	5.0	5.1	5.0	6.5	6.0
Eos%	%	6.6	5.0	9.7	6.0	10.7	8.0
Bas%	%	1.3	10.0	1.0	10.0	1.3	10.0
Neu#	10 ⁹ /L	1.68	0.50	4.36	0.70	11.77	2.00
Lym#	10 ⁹ /L	1.36	0.70	2.52	0.80	3.74	2.00
Mon#	10 ⁹ /L	0.14	0.30	0.41	0.50	1.22	1.10
Eos#	10 ⁹ /L	0.22	0.20	0.78	0.60	2.01	1.40
Bas#	10 ⁹ /L	0.04	0.40	0.08	0.80	0.24	1.80
RBC	10 ¹² /L	2.36	0.24	4.69	0.30	5.54	0.50
HGB	g/L	60	6	136	6	174	8
HCT	%	18.4	2.0	41.9	3.0	52.6	4.0
MCV	fL	77.9	5.0	89.4	5.0	95.0	6.0
MCH	pg	25.4	2.5	29.0	2.5	31.4	2.5
MCHC	g/L	327	30	324	30	331	30
RDW-CV	%	17.3	3.0	15.5	3.0	15.4	3.0
PLT	10 ⁹ /L	54	20	255	40	483	60
MPV	fL	8.5	3.0	8.2	3.0	8.3	3.0

No mean values and ranges are provided for research parameters-RDW-SD, PDW, PCT, P-LCR, and P-LCC. User may establish these based on their usage.



1	2	3	4	5
bg Трайност на отворен флакон: 14 дни	Параметър	Единица	Средно	Обхват
cs Stabilita otevřené ampulky: 14 dní	Parametr	Jednotka	Průměr	Rozsah
da Åben rør stabilitet: 14 dage	Parameter	Enhed	Middelværdi	Område
de Stabilität geöffneter Flaschen: 14 Tage	Parameter	Einheit	Mittelwert	Bereich
el Διατήρηση δείγματος μετά την αποσφράγιση: 14 ημέρα	Παράμετρος	Μονάδα	Μέση	Αναμενόμενο εύρος α
es Estabilidad de la cápsula abierta: 14 días	Parámetro	Unidad	Media	Intervalo
et Stabiilsus avatud vialis: 14 päevade arv	Parameeter	Ühik	Keskmine	Vahemik
fr Stabilité en flacon ouvert: 14 jour	Paramètre	Unité	Moyenne	Intervalle
hr Stabilnost otvorene bočice: 14 dana	Parametar	Jedinica	Sredstvo	Raspon
it Stabilità della fiala aperta: 14 giorni	Parametro	Unità	Media	Intervallo
lt Stabilumas atidarius buteliuka: 14 dienų	Parametras	Vienetas	Vidurkis	Intervalas
lv Atvērt flakonu stabilitāti: 14 dienas	Parametrs	Vienība	Nozīmē	Diapazons
no Åpen rør/glass stabilitet: 14 dager	Parameter	Enhhet	Middelverdi	Område
pl Trwałość po otwarciu fiolek: 14 Liczba dni	Kryterium	Jednostka	Średnia	Zakres
pt Estabilidade após abertura do frasco: 14 dias	Parâmetro	Unidade	Média	Intervalo
ro Termenul de valabilitate al fiolei desfăcută: 14 zile	Parametru	Unitate	Medie	Interval
ru Стабильность вскрытого флакона: 14 дней	Параметр	Единица измерения	Среднее	диапазон
sk Stabilita otvorenej injekčnej liekovky: 14 dní	Parameter	Jednotka	Priemer	Rozsah
sl Stabilnost odprte stekleničke: 14 dni	Parameter	Enota	Povprečno	Območje
sr Proizvod stabilan nakon otvaranja pakovanja: 14 dana	Parametar	Jedinica	Srednja vrednost	Opseg
sv Hållbarhet för öppnad flaska: 14 dagar	Parameter	Enhhet	Medelvärde	Intervall
tr Açık şişe stabilitesi gün: 14 gün	Parametre	Birim	Ortalama	Aralığı

31340-12
1504546



www.boule.com

Boule
Boule Medical AB
Domnarvsgatan 4
SE-163 53 Spånga, Sweden