

# Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22406-1K	2024-11-01
Low	22406-11	2024-11-01
Normal	22406-12	2024-11-01
High	22406-13	2024-11-05

Open vial stability 14 days <sup>1</sup>

MEDONIC	CONTROL		L
	LOT		22406-11
Medonic M-series/M32		Mean <sup>2</sup>	± Range <sup>3</sup>
RBC	10 <sup>12</sup> /L	2.21	0.12
MCV	fL	64.2	5.0
RDW	%	14.8	4.0
Hct	%	14.2	2.5
Plt	10 <sup>9</sup> /L	78	15
MPV	fL	11.1	1.5
WBC	10 <sup>9</sup> /L	3.4	0.3
Hgb	g/dL	5.3	0.3
MCH	pg	24.0	1.5
MCHC	g/dL	37.3	3.0
Lymph#	10 <sup>9</sup> /L	1.5	0.4
Mid#	10 <sup>9</sup> /L	0.2	0.3
Gran#	10 <sup>9</sup> /L	1.7	0.8
Lymph	%	44.5	7.0
Mid	%	5.4	5.0
Gran	%	50.1	7.0

CONTROL		N
LOT		22406-12
	Mean <sup>2</sup>	± Range <sup>3</sup>
	4.15	0.18
	80.6	5.0
	14.2	4.0
	33.4	3.0
	233	30
	11.6	1.5
	8.6	0.6
	11.9	0.4
	28.7	1.5
	35.6	3.0
	3.9	0.6
	0.6	0.5
	4.1	1.0
	45.0	5.0
	7.4	5.0
	47.6	5.0

CONTROL		H
LOT		22406-13
	Mean <sup>2</sup>	± Range <sup>3</sup>
	5.04	0.22
	92.4	5.0
	14.7	4.0
	46.6	3.5
	530	60
	11.7	1.5
	21.3	1.8
	16.2	0.5
	32.1	1.5
	34.8	3.0
	12.9	1.5
	1.2	1.0
	7.2	2.5
	60.5	5.0
	5.6	5.0
	33.9	5.0

1

bg	Трайност на отворен флакон 14 дни (bg)
hr	Stabilnost otvorene bočice 14 dani (hr)
cs	Stabilita otevřené ampulky 14 dní (cs)
da	Åben rør stabilitet 14 dage (da)
et	Stabiilsus avatud vialis 14 päevade (et)
fr	Stabilité en flacon ouvert 14 jour (fr)
de	Stabilität geöffneter Flaschen 14 tage (de)
el	Διάρθρωση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it	Stabilità della fiala aperta 14 giorni (it)
lv	Atvērt flakonu stabilitāti 14 dienas (lv)
lt	Stabilumas atidarius buteliuką 14 dienas (lt)
no	Åpen rør/glass stabilitet 14 dager (no)
pl	Trwałość po otwarciu fiołki 14 Liczba dni (pl)
pt	Estabilidade após abertura do frasco 14 dias (pt)
ro	Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru	Стабильность вскрытого флакона 14 дней (ru)
sr	Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk	Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl	Stabilita otvorenej injekčnej liekovky 14 dni (sl)
es	Estabilidad de la cápsula abierta 14 días (es)
sv	Hållbarhet för öppnad flaska 14 dagar (sv)
tr	Açık şişe stabilitesi 14 gün (tr)

2

Средно
Sredstvo
Průměr
Middelværdi
Keskmine
Moyenne
Mittelwert
Μέση
Media
Nozīmē
Vidurkis
Middelverdi
Średnia
Média
Medie
Среднее
Srednja vrednost
Priemer
Povprečno
Media
Medelvärde
Ortalama

3

Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervall
Aralıqı

# For Medonic M-series / M32 with barcode reader

## M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

## M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series  
 Lot: 22406-11  
 Prod.: Con-Diff Low  
 ID: 2240611+  
 Exp.: 1-Nov-2024

Instr.: M-series  
 Lot: 22406-12  
 Prod.: Con-Diff Norm  
 ID: 2240612+  
 Exp.: 1-Nov-2024

Instr.: M-series  
 Lot: 22406-13  
 Prod.: Con-Diff High  
 ID: 2240613+  
 Exp.: 5-Nov-2024

