

Boule Con-Diff

CONTROL

IVD



| Control | Lot | |
|-----------|----------|------------|
| Tri-Level | -- | -- |
| Low | 22406-41 | 2024-11-01 |
| Normal | 22406-42 | 2024-11-01 |
| High | 22406-43 | 2024-11-05 |

Open vial stability 14 days ¹

SWELAB

| Swelab Alfa/Alfa Plus | CONTROL | | L |
|-----------------------|---------------------|----------------------|----------|
| | LOT | 22406-41 | 22406-41 |
| | Mean ² | ± Range ³ | |
| RBC | 10 ¹² /L | 2.23 | 0.12 |
| MCV | fL | 66.4 | 5.0 |
| RDW | % | 14.5 | 4.0 |
| Hct | % | 14.8 | 2.5 |
| Plt | 10 ⁹ /L | 80 | 15 |
| MPV | fL | 12.3 | 1.5 |
| WBC | 10 ⁹ /L | 3.5 | 0.3 |
| Hgb | g/dL | 5.5 | 0.3 |
| MCH | pg | 24.7 | 1.5 |
| MCHC | g/dL | 37.2 | 3.0 |
| Lymph# | 10 ⁹ /L | 1.7 | 0.4 |
| Mid# | 10 ⁹ /L | 0.1 | 0.3 |
| Gran# | 10 ⁹ /L | 1.7 | 0.8 |
| Lymph | % | 48.5 | 7.0 |
| Mid | % | 4.0 | 5.0 |
| Gran | % | 47.5 | 7.0 |

| CONTROL | N | |
|---------|-------------------|----------------------|
| | LOT | 22406-42 |
| | Mean ² | ± Range ³ |
| | 4.15 | 0.18 |
| | 82.4 | 5.0 |
| | 12.9 | 4.0 |
| | 34.2 | 3.0 |
| | 239 | 30 |
| | 12.6 | 1.5 |
| | 8.6 | 0.6 |
| | 11.9 | 0.4 |
| | 28.7 | 1.5 |
| | 34.8 | 3.0 |
| | 4.2 | 0.6 |
| | 0.4 | 0.5 |
| | 4.0 | 1.0 |
| | 48.9 | 5.0 |
| | 4.8 | 5.0 |
| | 46.3 | 5.0 |

| CONTROL | H | |
|---------|-------------------|----------------------|
| | LOT | 22406-43 |
| | Mean ² | ± Range ³ |
| | 5.03 | 0.22 |
| | 94.2 | 5.0 |
| | 13.4 | 4.0 |
| | 47.4 | 3.5 |
| | 546 | 60 |
| | 12.9 | 1.5 |
| | 21.2 | 1.8 |
| | 16.1 | 0.5 |
| | 32.0 | 1.5 |
| | 34.0 | 3.0 |
| | 13.5 | 1.5 |
| | 0.7 | 1.0 |
| | 7.0 | 2.5 |
| | 63.7 | 5.0 |
| | 3.5 | 5.0 |
| | 32.8 | 5.0 |

1

| | |
|----|---|
| bg | Трайност на отворен флакон 14 дни (bg) |
| hr | Stabilnost otvorene bočice 14 dani (hr) |
| cs | Stabilita otevřené ampulky 14 dní (cs) |
| da | Åben rør stabilitet 14 dage (da) |
| et | Stabiilsus avatud vialis 14 päevade (et) |
| fr | Stabilité en flacon ouvert 14 jour (fr) |
| de | Stabilität geöffneter Flaschen 14 tage (de) |
| el | Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el) |
| it | Stabilità della fiala aperta 14 giorni (it) |
| lv | Atvērt flakonu stabilitāti 14 dienas (lv) |
| lt | Stabilumas atidarius buteliuką 14 dienas (lt) |
| no | Åpen rør/glass stabilitet 14 dager (no) |
| pl | Trwałość po otwarciu fiołki 14 Liczba dni (pl) |
| pt | Estabilidade após abertura do frasco 14 dias (pt) |
| ro | Termenul de valabilitate al fiolei desfăcute 14 zile (ro) |
| ru | Стабильность вскрытого флакона 14 дней (ru) |
| sr | Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr) |
| sk | Stabilita otvorenej injekčnej liekovky 14 dní (sk) |
| sl | Stabilita otvorenej injekčnej liekovky 14 dni (sl) |
| es | Estabilidad de la cápsula abierta 14 días (es) |
| sv | Hållbarhet för öppnad flaska 14 dagar (sv) |
| tr | Açık şişe stabilitesi 14 gün (tr) |

2

| | |
|--|------------------|
| | Средно |
| | Sredstvo |
| | Průměr |
| | Middelverdi |
| | Keskmine |
| | Moyenne |
| | Mittelwert |
| | Μέση |
| | Media |
| | Nozīmē |
| | Vidurkis |
| | Middelverdi |
| | Średnia |
| | Média |
| | Medie |
| | Среднее |
| | Srednja vrednost |
| | Priemer |
| | Povprečno |
| | Media |
| | Medelvärde |
| | Ortalama |

3

| | |
|--|---------------------|
| | Обхват |
| | Raspon |
| | Rozsah |
| | Område |
| | Vahemik |
| | Intervalle |
| | Bereich |
| | Αναμενόμενο εύρος α |
| | Intervallo |
| | Diapazons |
| | Intervalas |
| | Område |
| | Zakres |
| | Intervalo |
| | Interval |
| | диапазон |
| | Orseg |
| | Rozsah |
| | Območje |
| | Intervalo |
| | Intervall |
| | Aralığı |

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa
Lot: 22406-41
Prod.: Con-Diff Low
ID: 2240641+
Exp.: 1-Nov-2024

Instr.: Swelab Alfa
Lot: 22406-42
Prod.: Con-Diff Norm
ID: 2240642+
Exp.: 1-Nov-2024

Instr.: Swelab Alfa
Lot: 22406-43
Prod.: Con-Diff High
ID: 2240643+
Exp.: 5-Nov-2024

