

# Certificate of Analysis

<b>Product</b>	BOULE HYPOCHLORITE 4% CLEANER (100 mL)
<b>Product No</b>	1504385
<b>Lot No</b>	8012780724
<b>Manufacturing date</b>	2024-07-11
<b>Use before</b>	2026-01-31

<b>Test Hypochlorite</b>	<b>Specification</b>
Hypochlorite concentration	4 – <5 %
pH	11.5-13.5

Actual values meet specifications.

Date: 2024-08-01



*Ulrika Miefalk*  
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.  
All information in this document is property of Boule Medical AB.

