

LOT: 12408-52  
EXP: 2025-01-19  
Open Vial Stability: 16 days

EXIGO		Mean	± Range	Expected Range	
WBC	$10^9/L$	8.6	0.6	8.0	9.2
Hct	%	32.1	2.5	29.6	34.6
MCV	fL	77.9	5.0	72.9	82.9
RBC	$10^{12}/L$	4.12	0.18	3.94	4.30
Hgb	g/dL	11.4	0.4	11.0	11.8
MCH	pg	27.7	1.5	26.2	29.2
MCHC	g/dL	35.5	2.2	33.3	37.7
PLT	$10^9/L$	218	30	188	248

This assay sheet is also available on :

[www.boule.com](http://www.boule.com)

## For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo  
Lot: 12408-52  
Prod.: Vet-Con Norm  
ID: 1240852+  
Exp.: 19-Jan-2025

