

## CERTIFICATE OF ANALYSIS

**Product:** Swelab™AlfaDiluent

**Lot No:** 2403-004

**Product Part #:** 501-208


**Exp. Date:** 2027-02-28

**Manufacturing Date:** 03-29-2024

**Quality Testing Date:** 04-19-2024

Parameter	Specification	Result	Pass/Fail
pH (units):	6.50 – 6.90	6.76	Pass
Osmolality (mOsm/kg):	294 – 312	301	Pass
Conductivity (mS/cm):	16.98 – 17.93	17.26	Pass
Background (x10 <sup>3</sup> /μL):	0.0 – 8.0	2.3	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

  
\_\_\_\_\_  
Quality Assurance

04/19/24  
\_\_\_\_\_  
Date

