

Boule Con-Diff

CONTROL

IVD



Control	LOT	
Tri-Level	22407-6K	11/28/24
Low	22407-61	11/28/24
Normal	22407-62	12/02/24
High	22407-63	12/02/24

Open vial stability 14 days¹

SWELAB

CONTROL	L
LOT	22407-61

CONTROL	N
LOT	22407-62

CONTROL	H
LOT	22407-63

Swelab Alfa/Alfa Plus		Mean ²	± Range ³
RBC	10 ¹² /L	2.26	0.12
MCV	fL	67.6	5.0
RDW	%	14.4	4.0
Hct	%	15.3	2.5
Plt *	10 ⁹ /L	85	15
MPV	fL	12.1	1.5
WBC	10 ⁹ /L	3.2	0.3
Hgb	g/dl	5.6	0.3
MCH	pg	24.8	1.5
MCHC	g/dl	36.6	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.5	0.8
Lymph	%	49.2	7.0
Mid	%	4.1	5.0
Gran	%	46.7	7.0

		Mean ²	± Range ³
		4.22	0.18
		82.1	5.0
		13.2	4.0
		34.6	3.0
		240	30
		12.4	1.5
		8.7	0.6
		12.0	0.4
		28.4	1.5
		34.7	3.0
		4.4	0.6
		0.4	0.5
		3.9	1.0
		50.8	5.0
		4.3	5.0
		44.9	5.0

		Mean ²	± Range ³
		5.14	0.22
		92.9	5.0
		12.8	4.0
		47.8	3.5
		525	60
		12.7	1.5
		21.1	1.8
		16.1	0.5
		31.3	1.5
		33.7	3.0
		13.5	1.5
		0.8	1.0
		6.8	2.5
		63.9	5.0
		3.9	5.0
		32.2	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dny (cs)
 - da Åben r r stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 p evade (et)
 - fr Stabilit  en flacon ouvert 14 jour (fr)
 - de Stabilit t ge ffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilit  della fiala aperta 14 giorni (it)
 - lv Atv rt flakon u stabilit ti 14 dienas (lv)
 - lt Stabilitumas atidarius buteliuok  14 dienas (lt)
 - no Åpen r r/glass stabilitet 14 dager (no)
 - pl Trwało c po otwarciu fiolek 14 Liczba dni (pl)
 - pt Estabilidade ap s abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfacute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injek nej liekovky 14 dni (sk)
 - sl Stabilita otvorenej injek nej liekovky 14 dni (sl)
 - es Estabilidad de la c psula abierta 14 d as (es)
 - sv H llbarhet f r  ppnad flaska 14 dagar (sv)
 - tr A ık ŐiŐe stabilitesi 14 g n (tr)

- 2
- Средно
 - Sredstvo
 - Pr m r
 - Middelv rdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μ ση
 - Media
 - Noz m 
 - Vidurkis
 - Middelverdi
 - Srednia
 - M dia
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povpre no
 - Media
 - Medelv rde
 - Ortalama
- 3
- Обхват
 - Raspon
 - Rozsah
 - Omr de
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο ε ρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Omr de
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Obmo je
 - Intervalo
 - Intervall
 - Aral đı



<http://www.swelab.com/support/downloads>



1504054

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa
Lot: 22407-61
Prod.: Con-Diff Low
ID: 2240761+
Exp.: 28-Nov-2024

Instr.: Swelab Alfa
Lot: 22407-62
Prod.: Con-Diff Norm
ID: 2240762+
Exp.: 2-Dec-2024

Instr.: Swelab Alfa
Lot: 22407-63
Prod.: Con-Diff High
ID: 2240763+
Exp.: 2-Dec-2024

