

Certificate of Analysis


| | |
|---------------------------|---------------------------|
| Product | EXIGO EOS REAGENT (1.9 L) |
| Product No | 1504278 |
| Lot No | 2409-244 |
| Manufacturing Date | 2024-09-26 |
| Expiration Date | 2026-09-25 |

| Test | Specification |
|----------------------|----------------------|
| pH | 8.70 – 9.30 |
| Osmolality (mOsm/kg) | 125 - 145 |
| Conductivity (mS/cm) | 6.80 – 7.60 |
| EOS Performance * | Passes test |
| Blank EOS | ≤ 0.1 |

* Performance test using calibrated Blood cell counter at Boule Medical reference laboratory.

Actual values meet specifications.

Date: 2024-09-27



Pernilla Hammarstedt Rydberg
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

