

Certificate of Analysis

Product H50V DETERGENT, RFID
Product No 1440105 (5 L)
Lot No 24050417M4
Manufacturing Date 2024-05-22
Expiration Date 2026-05-21

| Test | Specification |
|------------|--|
| Appearance | Detergent should be transparent liquid without precipitation, particle or floccules |
| pH | 9.0 – 10.5 |
| Blank | WBC \leq 0.2 \times 10 ⁹ /L RBC \leq 0.02 \times 10 ¹² /L HGB \leq 1g/L PLT \leq 10 \times 10 ⁹ /L |

Actual values meet specifications.

Date: 2024-10-17



Ulrika Miefalk
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

