

# Boule Con-5Diff A1

CONTROL

IVD

LOT



Tri-Level	72409-1K	2025/01/05
Low	72409-11	2025/01/05
Normal	72409-12	2025/01/05
High	72409-13	2025/01/05

Open vial stability: 14 days <sup>1</sup>

## Medonic M51 / Swelab Lumi

Parameter <sup>2</sup>	Unit <sup>3</sup>	Low		Normal		High	
		Mean <sup>4</sup>	± Range <sup>5</sup>	Mean <sup>4</sup>	± Range <sup>5</sup>	Mean <sup>4</sup>	± Range <sup>5</sup>
WBC	10 <sup>9</sup> /L	3.34	± 0.50	7.97	± 1.00	18.65	± 2.50
Neu%	%	48.8	16.0	54.9	8.0	61.7	8.0
Lym%	%	38.4	21.0	30.3	8.0	19.6	12.0
Mon%	%	5.2	5.0	4.8	5.0	5.4	6.0
Eos%	%	7.6	5.0	10.0	6.0	13.3	8.0
Bas%	%	0.9	10.0	0.7	10.0	0.8	10.0
Neu#	10 <sup>9</sup> /L	1.63	0.50	4.39	0.70	11.51	2.00
Lym#	10 <sup>9</sup> /L	1.28	0.70	2.41	0.80	3.66	2.00
Mon#	10 <sup>9</sup> /L	0.17	0.30	0.38	0.50	1.00	1.10
Eos#	10 <sup>9</sup> /L	0.25	0.20	0.79	0.60	2.49	1.40
Bas#	10 <sup>9</sup> /L	0.03	0.40	0.05	0.80	0.14	1.80
RBC	10 <sup>12</sup> /L	2.29	0.24	4.58	0.30	5.43	0.50
HGB	g/L	58	6	132	6	170	8
HCT	%	18.4	2.0	41.2	3.0	52.5	4.0
MCV	fL	80.2	5.0	90.1	5.0	96.8	6.0
MCH	pg	25.3	2.5	28.8	2.5	31.3	2.5
MCHC	g/L	316	30	320	30	324	30
RDW-CV	%	15.5	3.0	15.0	3.0	14.6	3.0
PLT	10 <sup>9</sup> /L	58	20	251	40	454	60
MPV	fL	8.8	3.0	8.5	3.0	8.4	3.0

No mean values and ranges are provided for research parameters-RDW-SD, PDW, PCT, P-LCR, and P-LCC. User may establish these based on their usage.



<p><b>1</b> bg Трайност на отворен флакон: 14 дни cs Stabilita otevřené ampulky: 14 dny da Åben rør stabilitet: 14 dage de Stabilität geöffneter Flaschen: 14 tage el Διατήρηση δείγματος μετά την αποσφράγιση: 14 ημέρα es Estabilidad de la cápsula abierta: 14 días et Stabiilsus avatud vialis: 14 päevade arv fr Stabilité en flacon ouvert: 14 jour hr Stabilnost otvorene bočice: 14 dani it Stabilità della fiala aperta: 14 giorni lt Stabilumas atidarius buteliukoje: 14 dienu lv Atvērt flakonu stabilitāte: 14 dienas no Åpen rør/glass stabilitet: 14 dager pl Trwałość po otwarzeniu fiołki: 14 Liczba dni pt Estabilidade após abertura do frasco: 14 dias ro Termenul de valabilitate al fiolei desfăcute: 14 zile ru Стабильность вскрытого флакона: 14 дней sk Stabilita otvorenej injekčnej liekovky: 14 dní sl Stabilitost odprte stekleničke: 14 dni sr Proizvod stabilan nakon otvaranja pakovanja: 14 dana sv Hållbarhet för öppnad flaska: 14 dagar tr Acik sise stabilitesi qün: 14 qün</p>	<p><b>2</b> Параметър Parametr Parameter Parameter Parámetro Parámetro Paramètre Parametar Parametro Parametras Parameters Parameter Kryterium Parâmetro Parametru Параметр Parameter Parameter Parametar Parameter Parametre</p>	<p><b>3</b> Единица Jednotka Enhed Einheit Μονάδα Unidad Ühik Unité Jedinica Unità Vienetas Vienība Enhet Jednostka Unidade Unitate Единица измерения Jednotka Enota Jedinica Enhet Birim</p>	<p><b>4</b> Средно Průměr Middelværdi Mittelwert Μέση Media Keskmine Moyenne Sredstvo Media Vidurkis Nozīmē Middelverdi Średnia Média Medie Среднее Priemer Povprečno Srednja vrednost Medelvärde Ortalama</p>	<p><b>5</b> Обхват Rozsah Område Bereich Αναμενόμενο εύρος α Intervalo Vahemik Intervalle Raspon Intervallo Intervalas Diapazons Område Zakres Intervalo Interval диапазон Rozsah Območje Opseg Intervall Aralığı</p>
---	---	---	--	---

1504546 31340-12



[www.boule.com](http://www.boule.com)

**Boule**  
Boule Medical AB  
Domnarvsgatan 4  
SE-105 53 Spånga, Sweden