


 Open vial stability 14 days <sup>1</sup>

Control	Lot	
Tri-Level	22409-6K	01/31/25
Low	22409-61	01/31/25
Normal	22409-62	02/04/25
High	22409-63	02/04/25

**MEDONIC**

Medonic M-series/M32

		CONTROL LOT	L 22409-61
		Mean <sup>2</sup>	± Range <sup>3</sup>
RBC	10 <sup>12</sup> /L	2.30	0.12
MCV	fL	65.7	5.0
RDW	%	15.4	4.0
Hct	%	15.1	2.5
Plt	10 <sup>9</sup> /L	84	15
MPV	fL	11.1	1.5
WBC	10 <sup>9</sup> /L	3.4	0.3
Hgb	g/dL	5.6	0.3
MCH	pg	24.3	1.5
MCHC	g/dL	37.1	3.0
Lymph#	10 <sup>9</sup> /L	1.6	0.4
Mid#	10 <sup>9</sup> /L	0.2	0.3
Gran#	10 <sup>9</sup> /L	1.6	0.8
Lymph	%	45.7	7.0
Mid	%	5.8	5.0
Gran	%	48.5	7.0

	CONTROL LOT	N 22409-62	
		Mean <sup>2</sup>	± Range <sup>3</sup>
		4.14	0.18
		76.6	5.0
		14.4	4.0
		31.7	3.0
		235	30
		11.4	1.5
		8.6	0.6
		11.2	0.4
		27.1	1.5
		35.3	3.0
		4.0	0.6
		0.6	0.5
		4.0	1.0
		47.0	5.0
		6.8	5.0
		46.2	5.0

	CONTROL LOT	H 22409-63	
		Mean <sup>2</sup>	± Range <sup>3</sup>
		5.04	0.22
		94.4	5.0
		14.0	4.0
		47.6	3.5
		505	60
		11.9	1.5
		21.4	1.8
		16.3	0.5
		32.3	1.5
		34.2	3.0
		13.1	1.5
		1.3	1.0
		7.0	2.5
		61.3	5.0
		5.9	5.0
		32.8	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
  - hr Stabilnost otvorene bočice 14 dani (hr)
  - cs Stabilita otevřené ampulky 14 dny (cs)
  - da Åben rr stabilitet 14 dage (da)
  - et Stabiilsus avatud viaalis 14 pevade (et)
  - fr Stabilit en flacon ouvert 14 jour (fr)
  - de Stabilitt geffneter Flaschen 14 tage (de)
  - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
  - it Stabilit della fiala aperta 14 giorni (it)
  - lv Atvrt flakonu stabilitti 14 dienas (lv)
  - lt Stabilumas atidarius buteliuk 14 dienas (lt)
  - no Åpen rr/glass stabilitet 14 dager (no)
  - pl Trwaość po otwarciu fiołki 14 Liczba dni (pl)
  - pt Estabilidade aps abertura do frasco 14 dias (pt)
  - ro Termenul de valabilitate al fiolei desfcute 14 zile (ro)
  - ru Стабильность вскрытого флакона 14 дней (ru)
  - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
  - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
  - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
  - es Estabilidad de la cpsula abierta 14 dias (es)
  - sv Hllbarhet fr ppnad flaska 14 dagar (sv)
  - tr Aık ŐiŐe stabilitesi 14 gn (tr)

- 2
- Средно
  - Sredstvo
  - Prmr
  - Middelvrdi
  - Keskmine
  - Moyenne
  - Mittelwert
  - Μση
  - Media
  - Nozm
  - Vidurkis
  - Middelverdi
  - Srednia
  - Mdia
  - Medie
  - Среднее
  - Srednja vrednost
  - Priemer
  - Povprečno
  - Media
  - Medelvrde
  - Ortalama

- 3
- Обхват
  - Raspon
  - Rozsah
  - Omrde
  - Vahemik
  - Intervalle
  - Bereich
  - Αναμενμενο ερος α
  - Intervallo
  - Diapazons
  - Intervalas
  - Omrde
  - Zakres
  - Intervalo
  - Interval
  - диапазон
  - Opseg
  - Rozsah
  - Območje
  - Intervalo
  - Intervall
  - Aralgi



[http://www.medonic.se/  
support/downloads](http://www.medonic.se/support/downloads)



1504054

## For Medonic M-series / M32 with barcode reader


**M-Series:**

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

**M32:**

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series	Instr.: M-series	Instr.: M-series
Lot: 22409-61	Lot: 22409-62	Lot: 22409-63
Prod.: Con-Diff Low	Prod.: Con-Diff Norm	Prod.: Con-Diff High
ID: 2240961+	ID: 2240962+	ID: 2240963+
Exp.: 31-Jan-2025	Exp.: 4-Feb-2025	Exp.: 4-Feb-2025

<p>1:  Bou l e C B A V 2</p>	<p> Bou l e C B A V 2</p>	<p> Bou l e C B A V 2</p>
<p>2:  2240961+</p>	<p> 2240962+</p>	<p> 2240963+</p>
<p>3:  25013100001153030000</p>	<p> 25020400002153030000</p>	<p> 25020400003153030000</p>
<p>4:  10004031200000000000</p>	<p> 10004031400000000000</p>	<p> 10004030400000000000</p>
<p>5:  23012657500560315125</p>	<p> 41418766501120431730</p>	<p> 50422944501630547635</p>
<p>6:  08415111152431537130</p>	<p> 23530114152711535330</p>	<p> 50560119153231534230</p>
<p>7:  03403016040020301608</p>	<p> 08606040060060504010</p>	<p> 21418131150131007025</p>
<p>8:  15440457700585048570</p>	<p> 14440470500685046250</p>	<p> 14040613500595032850</p>
<p>9:  11119399140929125486</p>	<p> 11117575823428792685</p>	<p> 11117221422838222495</p>