


Boule Cal

CAL

IVD

Calibrator	Lot	
	22409-24	2024-12-15

Open vial stability 5 days

Medonic

Parameter	RBC $10^{12}/L$	MCV fL	Plt $10^9/L$	MPV fL	WBC $10^9/L$	Hgb g/dL
Medonic M-series / M32	4.13	75.4	217	11.0	8.4	11.2
Range \pm	0.10	2.0	15	1.5	0.3	0.2



<http://www.medonic.se/>
www.swelab.com/



Produced for Boule Medical AB
by Clinical Diagnostic Solutions



1800 NW 65th Avenue
Plantation, FL 33313

For Medonic M-series with barcode reader

1. In the Main Menu go to “Q/C Menu” “Enter Con/Cal” “Control/Calibrator Assay Value Input”

2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: M-series

Lot: 22409-24

Prod.: Cal

ID: 2240924+

Exp.: 15-Dec-2024



This bar-code sheet is also available at :

www.medonic.se
www.swelab.com
www.quintus5part.com