

# Boule Con-Diff

CONTROL

IVD



Control	LOT	
Tri-Level	22409-6K	01/31/25
Low	22409-61	01/31/25
Normal	22409-62	02/04/25
High	22409-63	02/04/25

Open vial stability 14 days<sup>1</sup>

## SWELAB

CONTROL	L
LOT	22409-61

CONTROL	N
LOT	22409-62

CONTROL	H
LOT	22409-63

### Swelab Alfa/Alfa Plus

		Mean <sup>2</sup>	± Range <sup>3</sup>
RBC	10 <sup>12</sup> /L	2.32	0.12
MCV	fL	67.5	5.0
RDW	%	15.1	4.0
Hct	%	15.7	2.5
Plt *	10 <sup>9</sup> /L	86	15
MPV	fL	12.2	1.5
WBC	10 <sup>9</sup> /L	3.4	0.3
Hgb	g/dl	5.7	0.3
MCH	pg	24.6	1.5
MCHC	g/dl	36.3	3.0
Lymph#	10 <sup>9</sup> /L	1.7	0.4
Mid#	10 <sup>9</sup> /L	0.1	0.3
Gran#	10 <sup>9</sup> /L	1.6	0.8
Lymph	%	49.0	7.0
Mid	%	4.2	5.0
Gran	%	46.8	7.0

	Mean <sup>2</sup>	± Range <sup>3</sup>
	4.13	0.18
	78.3	5.0
	13.4	4.0
	32.3	3.0
	241	30
	12.4	1.5
	8.7	0.6
	11.2	0.4
	27.1	1.5
	34.7	3.0
	4.4	0.6
	0.4	0.5
	3.9	1.0
	50.7	5.0
	4.3	5.0
	45.0	5.0

	Mean <sup>2</sup>	± Range <sup>3</sup>
	5.05	0.22
	96.0	5.0
	12.9	4.0
	48.5	3.5
	526	60
	12.9	1.5
	21.2	1.8
	16.3	0.5
	32.3	1.5
	33.6	3.0
	13.5	1.5
	0.8	1.0
	6.9	2.5
	63.9	5.0
	3.6	5.0
	32.5	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
  - hr Stabilitet otvorene bočice 14 dni (hr)
  - cs Stabilita otevřené ampulky 14 dní (cs)
  - da Åben r r stabilitet 14 dage (da)
  - et Stabiilsus avatud vialis 14 p evade (et)
  - fr Stabilit  en flacon ouvert 14 jour (fr)
  - de Stabilit t ge ffneter Flaschen 14 tage (de)
  - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημε ρα (el)
  - it Stabilit  della fiala aperta 14 giorni (it)
  - lv Atv rt flakonu stabilit ti 14 dienas (lv)
  - lt Stabilumas atidarius buteliuk  14 dienas (lt)
  - no Åpen r r/glass stabilitet 14 dager (no)
  - pl Trwa o c po otwarciu fioleki 14 Liczba dni (pl)
  - pt Estabilidade ap s abertura do frasco 14 dias (pt)
  - ro Termenul de valabilitate al fiolei desfacute 14 zile (ro)
  - ru Стабильность вскрытого флакона 14 дней (ru)
  - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
  - sk Stabilita otvorenej injek nej liekovky 14 dni (sk)
  - sl Stabilita otvorenej injek nej liekovky 14 dni (sl)
  - es Estabilidad de la c psula abierta 14 dias (es)
  - sv H llbarhet f r  ppnad flaska 14 dagar (sv)
  - tr A ık ŐiŐe stabilitesi 14 g n (tr)

- 2
- Средно
  - Sredstvo
  - Pr m r
  - Middelv rdi
  - Keskmine
  - Moyenne
  - Mittelwert
  - Μ ση
  - Media
  - Noz m 
  - Vidurkis
  - Middelverdi
  - Srednia
  - M dia
  - Medie
  - Среднее
  - Srednja vrednost
  - Priemer
  - Povpre no
  - Media
  - Medelv rde
  - Ortalama

- 3
- Обхват
  - Raspon
  - Rozsah
  - Omr de
  - Vahemik
  - Intervalle
  - Bereich
  - Αναμενόμενο ε ρος α
  - Intervallo
  - Diapazons
  - Intervalas
  - Omr de
  - Zakres
  - Intervalo
  - Interval
  - диапазон
  - Opseg
  - Rozsah
  - Obmo je
  - Intervalo
  - Intervall
  - Aralıđı



[http://www.swelab.com/  
support/downloads](http://www.swelab.com/support/downloads)



1504054

## For Swelab Alfa / Alfa Plus with barcode reader

### Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

### Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa  
Lot: 22409-61  
Prod.: Con-Diff Low  
ID: 2240961+  
Exp.: 31-Jan-2025

Instr.: Swelab Alfa  
Lot: 22409-62  
Prod.: Con-Diff Norm  
ID: 2240962+  
Exp.: 4-Feb-2025

Instr.: Swelab Alfa  
Lot: 22409-63  
Prod.: Con-Diff High  
ID: 2240963+  
Exp.: 4-Feb-2025

