

Boule Con-5Diff G2

CONTROL

IVD



| Control | Lot |
|-----------|----------|
| Tri-Level | 62410-1K |
| Low | 62410-11 |
| Normal | 62410-12 |
| High | 62410-13 |

| |
|------------|
| 2025-03-04 |
| 2025-03-04 |
| 2025-03-04 |
| 2025-03-04 |

Open vial stability:

14 days¹

Quintus

| WBC | 10 ⁹ /L |
|-------|---------------------|
| NEU# | 10 ⁹ /L |
| LYM# | 10 ⁹ /L |
| MON# | 10 ⁹ /L |
| EO# | 10 ⁹ /L |
| BAS# | 10 ⁹ /L |
| NEU | % |
| LYM | % |
| MON | % |
| EO | % |
| BAS | % |
| RBC | 10 ¹² /L |
| Hgb | g/dl |
| Hct | % |
| MCV | fL |
| MCH | pg |
| MCHC | g/dl |
| RDWcv | % |
| Plt | 10 ⁹ /L |
| MPV | fL |

| LOT L | 62410-11 |
|-------------------|----------------------|
| Mean ² | ± Range ³ |
| 4.1 | 0.5 |
| 2.3 | 0.5 |
| 1.3 | 0.4 |
| 0.3 | 0.3 |
| 0.2 | 0.3 |
| 0.0 | 0.2 |
| 55.4 | 9.0 |
| 32.6 | 9.9 |
| 7.0 | 8.7 |
| 4.8 | 4.5 |
| 0.2 | 1.0 |
| 2.22 | 0.15 |
| 6.1 | 0.5 |
| 17.1 | 2.5 |
| 77.2 | 5.5 |
| 27.5 | 2.8 |
| 35.7 | 3.6 |
| 17.4 | 5.0 |
| 80 | 25 |
| 9.8 | 3.0 |

| LOT N | 62410-12 |
|-------------------|----------------------|
| Mean ² | ± Range ³ |
| 9.0 | 0.8 |
| 4.3 | 1.3 |
| 3.6 | 0.8 |
| 0.7 | 0.7 |
| 0.4 | 0.3 |
| 0.0 | 0.2 |
| 47.4 | 8.0 |
| 40.1 | 9.9 |
| 7.5 | 8.5 |
| 4.9 | 3.0 |
| 0.1 | 1.0 |
| 4.10 | 0.18 |
| 13.0 | 0.6 |
| 35.3 | 3.0 |
| 86.2 | 5.5 |
| 31.7 | 2.8 |
| 36.8 | 3.0 |
| 14.3 | 5.0 |
| 245 | 40 |
| 10.0 | 3.0 |

| LOT H | 62410-13 |
|-------------------|----------------------|
| Mean ² | ± Range ³ |
| 20.7 | 2.0 |
| 8.0 | 3.4 |
| 10.6 | 2.7 |
| 0.8 | 1.9 |
| 1.3 | 0.6 |
| 0.0 | 0.2 |
| 38.5 | 8.0 |
| 50.8 | 9.9 |
| 4.1 | 8.3 |
| 6.5 | 3.6 |
| 0.1 | 1.0 |
| 4.92 | 0.30 |
| 17.7 | 0.7 |
| 45.7 | 3.5 |
| 92.8 | 5.5 |
| 36.0 | 2.8 |
| 38.7 | 3.0 |
| 15.0 | 5.0 |
| 518 | 65 |
| 10.4 | 3.0 |

| | CONTROL | N |
|---|----------|----------|
| i | Target X | Target Y |
| | 91 | 118 |

| | | |
|----|---|--|
| bg | 1 | Трайност на отворен флакон |
| hr | | Stabilnost otvorene bočice |
| cs | | Stabilita otevřené ampulky |
| da | | Åben r r stabilitet |
| et | | Stabiilsus avatud vialais |
| fr | | Stabilit  en flacon ouvert |
| de | | Stabilit t ge ffneter Flaschen |
| el | | Διατ ρηση δειγματος μετ  την αποσφραγιση |
| it | | Stabilit  della fiala aperta |
| lv | | Atv rt flakonu stabilit ti |
| lt | | Stabilumas atidarius buteliuk  |
| no | | Åpen r r/glass stabilitet |
| pl | | Trwa ość po otworzeniu fiolki |
| pt | | Estabilidade ap s abertura do frasco |
| ro | | Termenul de valabilitate al fiolei desf cut  |
| ru | | Стабильность вскрытого флакона |
| sr | | Proizvod stabilan nakon otvaranja pakovanja |
| sk | | Stabilita otvorenej injekčnej liekovky |
| sl | | Stabilnost odprte stekleni ce |
| es | | Estabilidad de la c psula abierta |
| sv | | H llbarhet f r  ppnad flaska |
| tr | | Açık ŐiŐe stabilitesi |

| | | | | |
|----------------|---|------------------|---|---------------------|
| 14  ni | 2 | Средно | 3 | Обхват |
| 14 dani | | Sredstvo | | Raspon |
| 14 dny | | Pr m r | | Rozsah |
| 14 dage | | Middelv rdi | | Omr de |
| 14 p evade arv | | Keskmine | | Vahemik |
| 14 jour | | Moyenne | | Intervalle |
| 14 tage | | Mittelwert | | Bereich |
| 14  m ra | | M soj | | Αναμενόμενο ε ρος α |
| 14 giorni | | Media | | Intervallo |
| 14 dienas | | Noz m  | | Diapazons |
| 14 d. | | Vidurkis | | Intervalas |
| 14 dager | | Middelverdi | | Omr de |
| 14 Liczba dni | | Őrednia | | Zakres |
| 14 dias | | M dia | | Intervalo |
| 14 zile | | Medie | | Interval |
| 14  ней | | Среднее | | диапазон |
| 14 dana | | Srednja vrednost | | Opseg |
| 14 d n  | | Priemer | | Rozsah |
| 14 dni | | Povpre no | | Obmo je |
| 14 dias | | Media | | Intervalo |
| 14 dagar | | Medelv rde | | Intervall |
| 14 g n | | Ortlama | | Aralıęı |



www.boule.com

For Quintus with barcode reader

1. In the Main Menu go to "Calibration QC" / "Enter QC reference".
2. Highlight the first Barcode field, scan in each calibrator/control barcode of the column in sequence (1-9) and click "OK".
3. Repeat steps 1 and 2 for each lot to be entered.
- 4.

Instr.: Quintus
 Lot: 62410-11
 Prod.: Con-diff Low
 ID: 6241011+
 Exp.: 4-Mar-2025

Instr.: Quintus
 Lot: 62410-12
 Prod.: Con-diff Normal
 ID: 6241012+
 Exp.: 4-Mar-2025

Instr.: Quintus
 Lot: 62410-13
 Prod.: Con-diff High
 ID: 6241013+
 Exp.: 4-Mar-2025

