

Certificate of Analysis

Product MEDONIC M-SERIES DUALPACK (1 L Lyse / 5 L Diluent)
Product No 1504128
Lot No 2410-256
Manufacturing Date 2024-10-23
Expiration Date 2027-10-22

| Test | Specification Diluent | Specification Lyse |
|------------------------------|-----------------------|--------------------|
| pH | 6,90 – 7,10 | 6,80 – 7,20 |
| Osmolality (mOsm/kg) | 290 – 310 | 220 – 240 |
| Conductivity (mS/cm) | 17,50 – 18,00 | 14,00 – 15,00 |
| CBC Performance * | Passes test | Passes test |
| Particles (PLT/WBC channels) | ≤ 5 | ≤ 0,2 |

* Performance test using calibrated Blood cell counter at Boule Medical reference laboratory. All parameters except mid cells are tested. Accepted reference 95 % to 105 % of expected control blood sample values.

Actual values meet specifications.

Date: 2024-11-22



Caroline Sparf
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

