

# CERTIFICATE OF ANALYSIS

Product: Boule Cleaning Kit

Kit Lot No: 2048

Product Part #: 501-036

Kit Exp. Date: 2026-07-31


<b>Part Number: 503-085</b>			
<b>Boule Detergent Cleaner Lot # 2045</b>			
<b>Manufacturing Date: 2024-07-11</b>		<b>Quality Testing Date: 2024-07-15</b>	
Parameter	Specification	Result	Pass/Fail
pH (units):	7.60 – 8.40	8.23	Pass

<b>Part Number: 503-086</b>			
<b>Boule Hypochlorite (2%) Cleaner Lot # 2046</b>			
<b>Manufacturing Date: 2024-07-09</b>		<b>Quality Testing Date: 2024-07-10</b>	
Parameter	Specification	Result	Pass/Fail
pH (units):	11.80 – 13.00	12.97	Pass

<b>Part Number: 503-087</b>			
<b>Boule Enzymatic Cleaner Lot # 2047</b>			
<b>Manufacturing Date: 2024-07-10</b>		<b>Quality Testing Date: 2024-07-15</b>	
Parameter	Specification	Result	Pass/Fail
pH (units):	7.60 – 8.40	8.09	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

  
\_\_\_\_\_  
Quality and Regulatory

  
\_\_\_\_\_  
Date

