

## CERTIFICATE OF ANALYSIS

**Product:** Quintus 5-Part Diluent

**Lot No:** 2411-201

**Product Part #:** 501-292

**Exp. Date:** 2026-10-31

**Manufacturing Date:** 2024-11-12

**Quality Testing Date:** 2024-11-21

Parameter	Specification	Result	Pass/Fail
pH (units):	6.7 – 7.1	7.1	Pass
Osmolality (mOsm/kg):	320 – 350	329	Pass
Conductivity (mS/cm):	18.50 – 19.50	18.71	Pass
Background ( $\times 10^3/\mu\text{L}$ )	0 – 5	0	Pass
Functional Test	Pass	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

*Deywill*  
Quality and Regulatory

12/02/24  
Date

