

LOT: 12411-52
EXP: 2025-04-21
Open Vial Stability: 16 days

EXIGO		Mean	± Range	Expected Range	
WBC	$10^9/L$	8.4	0.6	7.8	9.0
Hct	%	31.6	2.5	29.1	34.1
MCV	fL	76.5	5.0	71.5	81.5
RBC	$10^{12}/L$	4.13	0.18	3.95	4.31
Hgb	g/dL	11.2	0.4	10.8	11.6
MCH	pg	27.1	1.5	25.6	28.6
MCHC	g/dL	35.4	2.2	33.2	37.6
PLT	$10^9/L$	213	30	183	243

This assay sheet is also available on :
www.boule.com

For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12411-52
Prod.: Vet-Con Norm
ID: 1241152+
Exp.: 21-Apr-2025

