

## CERTIFICATE OF ANALYSIS

**Product:** Medonic™ M-Series Diluent

**Lot No:** 2501-005

**Product Part #:** 501-211

**Exp. Date:** 2027-12-31

**Manufacturing Date:** 2025-01-28


**Quality Testing Date:** 2025-01-30

Parameter	Specification	Result	Pass/Fail
pH (units):	6.90 – 7.20	7.13	Pass
Osmolality (mOsm/kg):	290 – 310	293	Pass
Conductivity (mS/cm):	17.50 – 18.00	17.78	Pass
Background Test ( $\times 10^3/\mu\text{L}$ ):	0.0 – 5.0	2.0	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Quality and Regulatory



Date

