

CERTIFICATE OF ANALYSIS

Product: Quintus 5-Part Diluent

Lot No: 2501-204

Product Part #: 501-292

Exp. Date: 2026-12-31

Manufacturing Date: 2025-01-31

Quality Testing Date: 2025-02-05

Parameter	Specification	Result	Pass/Fail
pH (units):	6.7 – 7.1	7.1	Pass
Osmolality (mOsm/kg):	320 – 350	327	Pass
Conductivity (mS/cm):	18.50 – 19.50	18.66	Pass
Background (x10 ³ /μL)	0 – 8	1	Pass
Functional Test	Pass	Pass	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

McLazari
Quality and Regulatory

2/6/25
Date

