

## CERTIFICATE OF ANALYSIS

**Product:** Quintus 5-Part Lyse

**Lot No:** 2501-301

**Product Part #:** 501-291

**Exp. Date:** 2026-07-31


**Manufacturing Date:** 2025-01-23

**Quality Testing Date:** 2025-02-26

Parameter	Specification	Result	Pass/Fail
pH (units):	2.7 – 3.0	2.9	Pass
Osmolality (mOsm/kg):	35 – 44	35	Pass
Conductivity (µS/cm):	1,117 – 1,541	1,253	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

  
\_\_\_\_\_  
Quality Assurance

  
\_\_\_\_\_  
Date

