

# Certificate of Analysis

**Product** MEDONIC M51-L2 LYSE (500 mL)  
**Product No** 1504512  
**Lot No** 2024092701  
**Manufacturing Date** 2024-09-27  
**Expiration Date** 2026-09-26

| Test                                  | Specification |
|---------------------------------------|---------------|
| pH                                    | 4,0 – 4,4     |
| CBC Performance test*                 | Passes test   |
| Background WBC (x 10 <sup>9</sup> /L) | ≤ 0,2         |

\*Performance testing at a reference laboratory. WBC, HGB and WBC scattergrams tested and approved.

Actual values meet specification.

Date: 2025-03-05



*Elin Andersson*  
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.  
All information in this document is property of Boule Medical AB.

