

LOT: 12502-52
EXP: 2025-07-20
Open Vial Stability: 16 days

EXIGO		Mean	± Range	Expected Range	
WBC	$10^9/L$	8.5	0.6	7.9	9.1
Hct	%	31.4	2.5	28.9	33.9
MCV	fL	75.5	5.0	70.5	80.5
RBC	$10^{12}/L$	4.16	0.18	3.98	4.34
Hgb	g/dL	11.2	0.4	10.8	11.6
MCH	pg	26.9	1.5	25.4	28.4
MCHC	g/dL	35.7	2.2	33.5	37.9
PLT	$10^9/L$	219	30	189	249

This assay sheet is also available on :

www.boule.com

For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12502-52
Prod.: Vet-Con Norm
ID: 1250252+
Exp.: 20-Jul-2025

