
Certificate of Analysis

Product H50V Probe Cleaner
Product No 1440111 (100 mL)
Lot No 25020219M7
Manufacturing Date 2025-02-11
Expiration Date 2026-02-10

Test	Specification
Appearance	Transparent liquid without precipitation, particle or floccules
pH	12,00 - 14,00
Blank	WBC \leq 0.5 \times 10 ⁹ /L RBC \leq 0.1 \times 10 ¹² /L HGB \leq 5g/L PLT \leq 20 \times 10 ⁹ /L

Actual values meet specifications.

Date: 2025-03-13



Pernilla Hammarstedt Rydberg
Quality Assurance

This product is intended for Laboratory use only unless otherwise specifically stated.
All information in this document is property of Boule Medical AB.

