

MEDONIC M-SERIES HEMATOLOGY ANALYZER - S/N: _____

Laboratory Name: _____

Month: _____ **Year** _____

Reference: Medonic M-series User's Manual

DAILY QUALITY CONTROL AND MAINTENANCE (User Manual, section 6.1 and 8.1)

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Check Reagent Levels																															
Check Printer Paper																															
Background Count																															
Quality Control																															
Clean Probe(s) w/alcohol																															
Performed by:																															

MONTHLY MAINTENANCE (User Manual, Section 8.2)

Boule Cleaning Kit Reagents	Lot #	Bottle #	Date	Tech
Monthly Cleaning (Hypochlorite)		2		
Clot Prevention (Enzymatic)		1		

Background Limits

WBC	≤ 0.1
RBC	≤ 0.01
HGB	≤ 0.2
PLT	≤ 10

SIX MONTH MAINTENANCE (User Manual, Section 8.3)

Boule Cleaning Kit Procedure (part # CDS 501-036)	Date	Tech
Lot #		

CALIBRATION (User Manual, Section 7)

Boule Cal Calibrator (Part # CDS 501-606)	Date	Tech
Lot #		

CURRENT QC INFORMATION

Boule Con-Diff Tri Level Control (Part # CDS 501-_____)	Exp. Date	Tech
Current Lot #		
Parallel Lot #		

Comments: