



## Preventive Maintenance Protocol - BM8X0

Instrument type	Serial number	Software version	Date

Customer address	Customer ID	Company address

### Maintenance checklist (X = Checked/Done, N/A = Not Applicable, ! = Comment below)

1. Perform cleaning kit 1504111		13. Check HBG adjustment	
2. Backflush shear valve		14. Check battery condition	
3. Clean mixing beaker		15. Check RBC/WBC counting time	
4. Clean MPA		16. Check noise test	
5. Clean CT module		17. Check MPA module	
6. Clean/replace CT needle		18. Check CT module	
7. Check/replace tubings		19. Waste pump pressure (>500 mbar)	
8. Replace MPA O-ring		20. Check air pump (>250 mbar)	
9. Perform waste pump repair kit		21. Check electromagnetic valves	
10. Check reagent detectors		22. Verify background results	
11. Check blood detector		23. Verify control blood results	
12. Check mixing beaker detectors		24. Calibrate, re-run control (if needed)	

### Comments

--

Customer signature	Technician signature
Customer name	Technician name