

# Boule Con-5Diff A1

CONTROL

IVD

LOT



Tri-Level	72411-1K	2025/03/05
Low	72411-11	2025/03/05
Normal	72411-12	2025/03/05
High	72411-13	2025/03/05

Open vial stability: 14 days <sup>1</sup>

## Medonic M51 / Swelab Lumi

Parameter <sup>2</sup>	Unit <sup>3</sup>	Low		Normal		High	
		Mean <sup>4</sup>	± Range <sup>5</sup>	Mean <sup>4</sup>	± Range <sup>5</sup>	Mean <sup>4</sup>	± Range <sup>5</sup>
WBC	10 <sup>9</sup> /L	3.44	± 0.50	8.09	± 1.00	18.62	± 2.50
Neu%	%	48.8	16.0	55.1	8.0	63.0	8.0
Lym%	%	40.8	21.0	30.3	8.0	19.3	12.0
Mon%	%	3.6	5.0	4.4	5.0	5.9	6.0
Eos%	%	6.8	5.0	10.2	6.0	11.8	8.0
Bas%	%	0.9	10.0	0.8	10.0	0.6	10.0
Neu#	10 <sup>9</sup> /L	1.68	0.50	4.45	0.70	11.73	2.00
Lym#	10 <sup>9</sup> /L	1.40	0.70	2.45	0.80	3.59	2.00
Mon#	10 <sup>9</sup> /L	0.12	0.30	0.36	0.50	1.10	1.10
Eos#	10 <sup>9</sup> /L	0.23	0.20	0.83	0.60	2.20	1.40
Bas#	10 <sup>9</sup> /L	0.03	0.40	0.07	0.80	0.10	1.80
RBC	10 <sup>12</sup> /L	2.40	0.24	4.75	0.30	5.53	0.50
HGB	g/L	62	6	137	6	174	8
HCT	%	19.6	2.0	42.8	3.0	53.5	4.0
MCV	fL	81.6	5.0	90.2	5.0	96.8	6.0
MCH	pg	25.8	2.5	28.8	2.5	31.5	2.5
MCHC	g/L	317	30	320	30	325	30
RDW-CV	%	14.8	3.0	14.8	3.0	14.4	3.0
PLT	10 <sup>9</sup> /L	65	20	246	40	466	60
MPV	fL	8.6	3.0	7.9	3.0	8.5	3.0

No mean values and ranges are provided for research parameters-RDW-SD, PDW, PCT, P-LCR, and P-LCC. User may establish these based on their usage.



<p><b>1</b></p> <p>bg Трайност на отворен флакон: 14 дни          cs Stabilita otevřené ampulky: 14 dny          da Åben rør stabilitet: 14 dage          de Stabilität geöffneter Flaschen: 14 tage          el Διατήρηση δείγματος μετά την αποσφράγιση: 14 ημέρα          es Estabilidad de la cápsula abierta: 14 días          et Stabiilsus avatud vialis: 14 päevade arv          fr Stabilité en flacon ouvert: 14 jour          hr Stabilnost otvorene bočice: 14 dani          it Stabilità della fiala aperta: 14 giorni          lt Stabilumas atidarius buteliuko: 14 dienų          lv Atvērt flakonu stabilitāti: 14 dienas          no Åpen rør/ glass stabilitet: 14 dager          pl Trwałość po otworzeniu fiołki: 14 Liczba dni          pt Estabilidade após abertura do frasco: 14 dias          ro Termenul de valabilitate al fiolei desfăcută: 14 zile          ru Стабильность вскрытого флакона: 14 дней          sk Stabilita otvorenej injekčnej liekovky: 14 dní          sl Stabilitost odprte stekleničke: 14 dni          sr Proizvod stabilan nakon otvaranja pakovanja: 14 dana          sv Hållbarhet för öppnad flaska: 14 dagar          tr Acik sise stabilitesi qün: 14 qün</p>	<p><b>2</b></p> <p>Параметър          Parametr          Parameter          Parameter          Παράμετρος          Parámetro          Parameeter          Paramètre          Parametar          Parametro          Parametras          Parametrs          Parameter          Kryterium          Parâmetro          Parametru          Параметр          Parameter          Parameter          Parametar          Parameter          Parametre</p>	<p><b>3</b></p> <p>Единица          Jednotka          Enhed          Einheit          Μονάδα          Unidad          Ühik          Unité          Jedinica          Unità          Vienetas          Vienība          Enhet          Jednostka          Unidade          Unitate          Единица измерения          Jednotka          Enota          Jedinica          Enhet          Birim</p>	<p><b>4</b></p> <p>Средно          Průměr          Middelværdi          Mittelwert          Μέση          Media          Keskmine          Moyenne          Sredstvo          Media          Vidurkis          Nozīmē          Middelverdi          Średnia          Média          Médie          Среднее          Priemer          Povprečno          Srednja vrednost          Medelvärde          Ortalama</p>	<p><b>5</b></p> <p>Обхват          Rozsah          Område          Bereich          Αναμενόμενο εύρος α          Intervalo          Vahemik          Intervalle          Raspon          Intervallo          Intervalas          Diapazons          Område          Zakres          Intervalo          Intervall          диапазон          Rozsah          Območje          Opseg          Intervall          Aralıđı</p>
--	--	---	--	--

1504546 31340-12