

Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22503-1K	2025-08-04
Low	22503-11	2025-08-04
Normal	22503-12	2025-08-04
High	22503-13	2025-08-06

Open vial stability 14 days ¹

MEDONIC		CONTROL	L
Medonic M-series/M32		LOT	22503-11
		Mean ²	± Range ³
RBC	10 ¹² /L	2.28	0.12
MCV	fL	65.7	5.0
RDW	%	15.9	4.0
Hct	%	15.0	2.5
Plt	10 ⁹ /L	83	15
MPV	fL	11.1	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	24.1	1.5
MCHC	g/dL	36.7	3.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	43.1	7.0
Mid	%	6.0	5.0
Gran	%	50.9	7.0

CONTROL	N	
LOT	22503-12	
	Mean ²	± Range ³
	4.23	0.18
	82.8	5.0
	13.4	4.0
	35.0	3.0
	240	30
	11.4	1.5
	8.3	0.6
	12.3	0.4
	29.1	1.5
	35.1	3.0
	3.9	0.6
	0.6	0.5
	3.8	1.0
	46.8	5.0
	7.2	5.0
	46.0	5.0

CONTROL	H	
LOT	22503-13	
	Mean ²	± Range ³
	5.16	0.22
	93.1	5.0
	13.6	4.0
	48.0	3.5
	541	60
	11.6	1.5
	20.3	1.8
	16.5	0.5
	32.0	1.5
	34.4	3.0
	12.5	1.5
	1.2	1.0
	6.6	2.5
	61.6	5.0
	5.9	5.0
	32.5	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelværdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Interval
 - Intervall
 - Aralığı

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series	Instr.: M-series	Instr.: M-series
Lot: 22503-11	Lot: 22503-12	Lot: 22503-13
Prod.: Con-Diff Low	Prod.: Con-Diff Norm	Prod.: Con-Diff High
ID: 2250311+	ID: 2250312+	ID: 2250313+
Exp.: 4-Aug-2025	Exp.: 4-Aug-2025	Exp.: 6-Aug-2025

1:  Boul e CBAV2	 Boul e CBAV2	 Boul e CBAV2
2:  2250311+	 2250312+	 2250313+
3:  25080400001153030000	 25080400002153030000	 25080600003153030000
4:  10004030700000000000	 10004030200000000000	 10004030100000000000
5:  22812657500550315025	 42318828501230435030	 51622931501650548035
6:  08315111152411536730	 24030114152911535130	 54160116153201534430
7:  03403015040020301708	 08306039060060503810	 20318125150121006625
8:  15940431700605050970	 13440468500725046050	 13640616500595032550
9:  11110337629432786958	 11119729746994293323	 11118482841925630043

Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22503-1K	2025-08-04
Low	22503-11	2025-08-04
Normal	22503-12	2025-08-04
High	22503-13	2025-08-06

Open vial stability 14 days ¹

SWELAB

CONTROL		L	
LOT		22503-11	
Swelab Alfa/Alfa Plus		Mean ²	± Range ³
RBC	10 ¹² /L	2.31	0.12
MCV	fL	67.0	5.0
RDW	%	15.1	4.0
Hct	%	15.5	2.5
Plt	10 ⁹ /L	86	15
MPV	fL	11.8	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.7	0.3
MCH	pg	24.7	1.5
MCHC	g/dL	36.8	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	47.1	7.0
Mid	%	3.3	5.0
Gran	%	49.6	7.0

CONTROL		N	
LOT		22503-12	
		Mean ²	± Range ³
		4.29	0.18
		84.1	5.0
		12.4	4.0
		36.1	3.0
		249	30
		12.0	1.5
		8.3	0.6
		12.4	0.4
		28.9	1.5
		34.3	3.0
		4.3	0.6
		0.3	0.5
		3.7	1.0
		50.9	5.0
		4.1	5.0
		45.0	5.0

CONTROL		H	
LOT		22503-13	
		Mean ²	± Range ³
		5.21	0.22
		94.5	5.0
		12.6	4.0
		49.2	3.5
		569	60
		12.3	1.5
		20.2	1.8
		16.4	0.5
		31.5	1.5
		33.3	3.0
		13.1	1.5
		0.7	1.0
		6.4	2.5
		64.8	5.0
		3.5	5.0
		31.7	5.0

1

bg	Трайност на отворен флакон 14 дни (bg)
hr	Stabilnost otvorene bočice 14 dani (hr)
cs	Stabilita otevřené ampulky 14 dní (cs)
da	Åben rør stabilitet 14 dage (da)
et	Stabiilsus avatud viaalis 14 päevade (et)
fr	Stabilité en flacon ouvert 14 jour (fr)
de	Stabilität geöffneter Flaschen 14 tage (de)
el	Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it	Stabilità della fiala aperta 14 giorni (it)
lv	Atvērt flakonu stabilitāti 14 dienas (lv)
lt	Stabilumas atidarius buteliuką 14 dienas (lt)
no	Åpen rør/glass stabilitet 14 dager (no)
pl	Trwałość po otwarciu fiołki 14 Liczba dni (pl)
pt	Estabilidade após abertura do frasco 14 dias (pt)
ro	Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru	Стабильность вскрытого флакона 14 дней (ru)
sr	Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk	Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl	Stabilita otvorenej injekčneje liekovky 14 dni (sl)
es	Estabilidad de la cápsula abierta 14 días (es)
sv	Hållbarhet för öppnad flaska 14 dagar (sv)
tr	Açık şişe stabilitesi 14 gün (tr)

2

	Средно
	Sredstvo
	Průměr
	Middelværdi
	Keskmine
	Moyenne
	Mittelwert
	Μέση
	Media
	Nozīmē
	Vidurkis
	Middelverdi
	Średnia
	Média
	Medie
	Среднее
	Srednja vrednost
	Priemer
	Povprečno
	Media
	Medelvärde
	Ortalama

3

	Обхват
	Raspon
	Rozsah
	Område
	Vahemik
	Intervalle
	Bereich
	Αναμενόμενο εύρος α
	Intervallo
	Diapazons
	Intervalas
	Område
	Zakres
	Intervalo
	Interval
	диапазон
	Orseg
	Rozsah
	Območje
	Intervalo
	Intervall
	Aralığı

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa
Lot: 22503-11
Prod.: Con-Diff Low
ID: 2250311+
Exp.: 4-Aug-2025

Instr.: Swelab Alfa
Lot: 22503-12
Prod.: Con-Diff Norm
ID: 2250312+
Exp.: 4-Aug-2025

Instr.: Swelab Alfa
Lot: 22503-13
Prod.: Con-Diff High
ID: 2250313+
Exp.: 6-Aug-2025

