

## CERTIFICATE OF ANALYSIS

**Product:** Swelab™AlfaDiluent

**Lot No:** 2401-004

**Product Part #:** 501-208

**Exp. Date:** 2026-12-31

**Manufacturing Date:** 01-12-2024

**Quality Testing Date:** 01-17-2024

Parameter	Specification	Result	Pass/Fail
pH (units):	6.50 – 6.90	6.80	Pass
Osmolality (mOsm/kg):	294 – 312	299	Pass
Conductivity (mS/cm):	16.98 – 17.93	17.27	Pass
Background (x10 <sup>3</sup> /μL):	0.0 – 8.0	0.3	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.



Quality Assurance



Date

