

## CERTIFICATE OF ANALYSIS

**Product:** Swelab™ AlfaDiluent RFID

**Lot No:** 2410-001

**Product Part #:** 501-465


**Exp. Date:** 2027-09-30


**Manufacturing Date:** 10-24-2024

**Quality Testing Date:** 11-05-2024

Parameter	Specification	Result	Pass/Fail
pH (units):	6.50 – 6.90	6.90	Pass
Osmolality (mOsm/kg):	294 – 312	297	Pass
Conductivity (mS/cm):	16.98 – 17.93	17.20	Pass
Background (x10 <sup>3</sup> /μL):	0.0 – 8.0	3.5	Pass

The product has been manufactured in compliance with Quality Systems Regulations and has met all current release specifications.

  
\_\_\_\_\_  
Quality Assurance

  
\_\_\_\_\_  
Date

