

Boule Con-5Diff A1

CONTROL

IVD

LOT



Tri-Level	72503-1K	2025/07/05
Low	72503-11	2025/07/05
Normal	72503-12	2025/07/05
High	72503-13	2025/07/05

Open vial stability: 14 days ¹

Medonic M51 / Swelab Lumi

Parameter ²	Unit ³	Low		Normal		High	
		Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵
WBC	10 ⁹ /L	3.51	0.50	8.13	1.00	19.22	2.50
Neu%	%	50.7	16.0	55.7	8.0	63.3	8.0
Lym%	%	39.1	21.0	32.7	8.0	20.1	12.0
Mon%	%	4.2	5.0	3.7	5.0	6.5	6.0
Eos%	%	6.0	5.0	7.9	6.0	10.1	8.0
Bas%	%	2.3	10.0	1.7	10.0	2.5	10.0
Neu#	10 ⁹ /L	1.78	0.50	4.53	0.70	12.17	2.00
Lym#	10 ⁹ /L	1.37	0.70	2.66	0.80	3.87	2.00
Mon#	10 ⁹ /L	0.15	0.30	0.30	0.50	1.25	1.10
Eos#	10 ⁹ /L	0.21	0.20	0.64	0.60	1.93	1.40
Bas#	10 ⁹ /L	0.08	0.40	0.14	0.80	0.47	1.80
RBC	10 ¹² /L	2.29	0.24	4.74	0.30	5.51	0.50
HGB	g/L	61	6	135	6	174	8
HCT	%	19.2	2.0	43.2	3.0	54.1	4.0
MCV	fL	83.8	5.0	91.2	5.0	98.1	6.0
MCH	pg	26.6	2.5	28.5	2.5	31.6	2.5
MCHC	g/L	318	30	312	30	322	30
RDW-CV	%	15.2	3.0	13.8	3.0	13.5	3.0
PLT	10 ⁹ /L	66	20	246	40	491	60
MPV	fL	8.0	3.0	7.8	3.0	7.9	3.0

No mean values and ranges are provided for research parameters-RDW-SD, PDW, PCT, P-LCR, and P-LCC. User may establish these based on their usage.



1	2	3	4	5
bg Трайност на отворен флакон: 14 дни	Параметър	Единица	Средно	Обхват
cs Stabilita otevřené ampulky: 14 dní	Parametr	Jednotka	Průměr	Rozsah
da Åben rør stabilitet: 14 dage	Parameter	Enhed	Middelværdi	Område
de Stabilität geöffneter Flaschen: 14 tage	Parameter	Einheit	Mittelwert	Bereich
el Διατήρηση δείγματος μετά την αποσφράγιση: 14 ημέρα	Παράμετρος	Μονάδα	Μέση	Αναμενόμενο εύρος α
es Estabilidad de la cápsula abierta: 14 días	Parámetro	Unidad	Media	Intervalo
et Stabiilsus avatud vialis: 14 päevade arv	Parameeter	Ühik	Keskmine	Vahemik
fr Stabilité en flacon ouvert: 14 jour	Paramètre	Unité	Moyenne	Intervalle
hr Stabilnost otvorene bočice: 14 dani	Parametar	Jedinica	Sredstvo	Raspon
it Stabilità della fiala aperta: 14 giorni	Parametro	Unità	Media	Intervallo
lt Stabilumas atidarius buteliuko: 14 dienų	Parametras	Vienetas	Vidurkis	Intervalas
lv Atvērt flakonu stabilitāti: 14 dienas	Parametrs	Vienība	Nozīmē	Diapazons
no Åpen rør/lass stabilitet: 14 dager	Parameter	Enhhet	Middelverdi	Område
pl Trwałość po otwarciu fiolek: 14 Liczba dni	Kryterium	Jednostka	Srednia	Zakres
pt Estabilidade após abertura do frasco: 14 dias	Parâmetro	Unidade	Média	Intervalo
ro Termenul de valabilitate al fiolei destăcut: 14 zile	Parametru	Unitate	Medie	Interval
ru Стабильность вскрытого флакона: 14 дней	Параметр	Единица измерения	Среднее	диапазон
sk Stabilita otvorenej injekčnej liekovky: 14 dní	Parameter	Jednotka	Priemer	Rozsah
sl Stabilnost odprte stekleničke: 14 dni	Parameter	Enota	Povprečno	Območje
sr Proizvod stabilan nakon otvaranja pakovanja: 14 dana	Parametar	Jedinica	Srednja vrednost	Opseg
sv Hållbarhet för öppnad flaska: 14 dagar	Parameter	Enhhet	Medelvärde	Intervall
tr Açık sise stabilitesi gün: 14 gün	Parametre	Birim	Ortalama	Aralığı

1504546 31340-12



www.boule.com

Boule
Boule Medical AB
Domnarvsgatan 4
SE-105 53 Spånga, Sweden