

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	--	--
Low	22504-41	2025-08-28
Normal	22504-42	2025-08-28
High	22504-43	2025-09-01

MEDONIC

CONTROL		L	
LOT		22504-41	
Medonic M-series/M32		Mean ²	± Range ³
RBC	10 ¹² /L	2.35	0.12
MCV	fL	63.8	5.0
RDW	%	16.1	4.0
Hct	%	15.0	2.5
Plt	10 ⁹ /L	81	15
MPV	fL	11.1	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	23.4	1.5
MCHC	g/dL	36.7	3.0
Lymph#	10 ⁹ /L	1.4	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	42.9	7.0
Mid	%	6.4	5.0
Gran	%	50.7	7.0

CONTROL		N	
LOT		22504-42	
		Mean ²	± Range ³
		4.22	0.18
		83.5	5.0
		13.0	4.0
		35.2	3.0
		241	30
		11.1	1.5
		8.3	0.6
		12.3	0.4
		29.1	1.5
		34.9	3.0
		3.9	0.6
		0.5	0.5
		3.9	1.0
		46.9	5.0
		6.2	5.0
		46.9	5.0

CONTROL		H	
LOT		22504-43	
		Mean ²	± Range ³
		5.05	0.22
		89.2	5.0
		12.9	4.0
		45.0	3.5
		540	60
		11.2	1.5
		20.4	1.8
		15.5	0.5
		30.7	1.5
		34.4	3.0
		12.5	1.5
		1.2	1.0
		6.7	2.5
		61.7	5.0
		5.7	5.0
		32.6	5.0

1

bg	Трайност на отворен флакон 14 дни (bg)
hr	Stabilnost otvorene bočice 14 dni (hr)
cs	Stabilita otevřené ampulky 14 dní (cs)
da	Åben rør stabilitet 14 dage (da)
et	Stabiilsus avatud viaalis 14 päevade (et)
fr	Stabilité en flacon ouvert 14 jour (fr)
de	Stabilität geöffneter Flaschen 14 tage (de)
el	Διάρκηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it	Stabilità della fiala aperta 14 giorni (it)
lv	Atvērt flakonu stabilitāti 14 dienas (lv)
lt	Stabilumas atidarius buteliuką 14 dienas (lt)
no	Åpen rør/glass stabilitet 14 dager (no)
pl	Trwałość po otwarciu fiołki 14 Liczba dni (pl)
pt	Estabilidade após abertura do frasco 14 dias (pt)
ro	Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru	Стабильность вскрытого флакона 14 дней (ru)
sr	Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk	Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl	Stabilita otvorenej injekčnej liekovky 14 dni (sl)
es	Estabilidad de la cápsula abierta 14 días (es)
sv	Hållbarhet för öppnad flaska 14 dagar (sv)
tr	Açık şişe stabilitesi 14 gün (tr)

2

Sredno
Sredstvo
Průměr
Middelværdi
Keskmine
Moyenne
Mittelwert
Μέση
Media
Nozīmē
Vidurkis
Middelverdi
Średnia
Média
Medie
Среднее
Srednja vrednost
Priemer
Povprečno
Media
Medelvärde
Ortalama

3

Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervalo
Intervall
Aralığı

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series
Lot: 22504-41
Prod.: Con-Diff Low
ID: 2250441+
Exp.: 28-Aug-2025

Instr.: M-series
Lot: 22504-42
Prod.: Con-Diff Norm
ID: 2250442+
Exp.: 28-Aug-2025

Instr.: M-series
Lot: 22504-43
Prod.: Con-Diff High
ID: 2250443+
Exp.: 1-Sep-2025



3.

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	--	--
Low	22504-41	2025-08-28
Normal	22504-42	2025-08-28
High	22504-43	2025-09-01

SWELAB

CONTROL		L	
LOT		22504-41	
Swelab Alfa/Alfa Plus		Mean ²	± Range ³
RBC	10 ¹² /L	2.36	0.12
MCV	fL	65.3	5.0
RDW	%	15.4	4.0
Hct	%	15.4	2.5
Plt	10 ⁹ /L	82	15
MPV	fL	11.8	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.6	0.3
MCH	pg	23.7	1.5
MCHC	g/dL	36.4	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.6	0.8
Lymph	%	47.5	7.0
Mid	%	3.8	5.0
Gran	%	48.7	7.0

CONTROL		N	
LOT		22504-42	
		Mean ²	± Range ³
		4.21	0.18
		84.7	5.0
		11.9	4.0
		35.7	3.0
		250	30
		11.7	1.5
		8.4	0.6
		12.3	0.4
		29.2	1.5
		34.5	3.0
		4.2	0.6
		0.3	0.5
		3.9	1.0
		50.1	5.0
		3.7	5.0
		46.2	5.0

CONTROL		H	
LOT		22504-43	
		Mean ²	± Range ³
		5.03	0.22
		90.3	5.0
		11.7	4.0
		45.4	3.5
		556	60
		12.0	1.5
		20.4	1.8
		15.5	0.5
		30.8	1.5
		34.1	3.0
		13.2	1.5
		0.7	1.0
		6.5	2.5
		64.4	5.0
		3.6	5.0
		32.0	5.0

1

bg	Трайност на отворен флакон 14 дни (bg)
hr	Stabilitet otvorene bočice 14 dni (hr)
cs	Stabilita otevřené ampulky 14 dní (cs)
da	Åben rør stabilitet 14 dage (da)
et	Stabiilsus avatud vialis 14 päevade (et)
fr	Stabilité en flacon ouvert 14 jour (fr)
de	Stabilität geöffneter Flaschen 14 tage (de)
el	Διάρθρωση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it	Stabilità della fiala aperta 14 giorni (it)
lv	Atvērt flakonu stabilitāti 14 dienas (lv)
lt	Stabilumas atidarius buteliuką 14 dienas (lt)
no	Åpen rør/glass stabilitet 14 dager (no)
pl	Trwałość po otwarciu fiolek 14 dni (pl)
pt	Estabilidade após abertura do frasco 14 dias (pt)
ro	Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru	Стабильность вскрытого флакона 14 дней (ru)
sr	Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk	Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl	Stabilita otvorenej injekčneje liekovky 14 dni (sl)
es	Estabilidad de la cápsula abierta 14 días (es)
sv	Hållbarhet för öppnad flaska 14 dagar (sv)
tr	Açık şişe stabilitesi 14 gün (tr)

2

Sredno
Sredstvo
Průměr
Middelverdi
Keskmine
Moyenne
Mittelwert
Μέση
Media
Nožmē
Vidurkis
Middelverdi
Średnia
Média
Medie
Среднее
Srednja vrednost
Priemer
Povprečno
Media
Medelvärde
Ortalama

3

Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervalo
Intervall
Aralığı

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa

Lot: 22504-41

Prod.: Con-Diff Low

ID: 2250441+

Exp.: 28-Aug-2025

Instr.: Swelab Alfa

Lot: 22504-42

Prod.: Con-Diff Norm

ID: 2250442+

Exp.: 28-Aug-2025

Instr.: Swelab Alfa

Lot: 22504-43

Prod.: Con-Diff High

ID: 2250443+

Exp.: 1-Sep-2025

1:



B o u l e C B A V 3



B o u l e C B A V 3



B o u l e C B A V 3

2:



2250441+



2250442+



2250443+

3:



25082800001153030000



25082800002153030000



25090100003153030000

4:



20014018000000000000



20014018000000000000



20014018000000000000

5:



23612653500560315425



42118847501230435730



50322903501550545435

6:



08215118152371536430



25030117152921534530



55660120153081534130

7:



03303016040010301608



08406042060030503910



20418132150071006525

8:



15440475700385048770



11940501500375046250



11740644500365032050

9:



11110814890119833409



11110298637229842286



11110681393851889329