



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22504-6K	08/28/25
Low	22504-61	08/28/25
Normal	22504-62	08/28/25
High	22504-63	09/01/25

MEDONIC

Medonic M-series/M32

		CONTROL	L
		LOT	22504-61
		Mean ²	± Range ³
RBC	10 ¹² /L	2.35	0.12
MCV	fL	63.8	5.0
RDW	%	16.1	4.0
Hct	%	15.0	2.5
Plt	10 ⁹ /L	81	15
MPV	fL	11.1	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	23.4	1.5
MCHC	g/dL	36.7	3.0
Lymph#	10 ⁹ /L	1.4	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	42.9	7.0
Mid	%	6.4	5.0
Gran	%	50.7	7.0

		CONTROL	N
		LOT	22504-62
		Mean ²	± Range ³
		4.22	0.18
		83.5	5.0
		13.0	4.0
		35.2	3.0
		241	30
		11.1	1.5
		8.3	0.6
		12.3	0.4
		29.1	1.5
		34.9	3.0
		3.9	0.6
		0.5	0.5
		3.9	1.0
		46.9	5.0
		6.2	5.0
		46.9	5.0

		CONTROL	H
		LOT	22504-63
		Mean ²	± Range ³
		5.05	0.22
		89.2	5.0
		12.9	4.0
		45.0	3.5
		540	60
		11.2	1.5
		20.4	1.8
		15.5	0.5
		30.7	1.5
		34.4	3.0
		12.5	1.5
		1.2	1.0
		6.7	2.5
		61.7	5.0
		5.7	5.0
		32.6	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
- hr Stabilität otvorene bočice 14 dni (hr)
- cs Stabilita otevřené ampulky 14 dní (cs)
- da Åben rør stabilitet 14 dage (da)
- et Stabiilsus avatud viaalis 14 päevade (et)
- fr Stabilité en flacon ouvert 14 jour (fr)
- de Stabilität geöffneter Flaschen 14 tage (de)
- el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
- it Stabilità della fiala aperta 14 giorni (it)
- lv Atvērt flakonu stabilitāti 14 dienas (lv)
- lt Stabilumas atidarius buteliuką 14 dienas (lt)
- no Åpen rør/glass stabilitet 14 dager (no)
- pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
- pt Estabilidade após abertura do frasco 14 dias (pt)
- ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
- ru Стабильность вскрытого флакона 14 дней (ru)
- sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
- sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
- sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
- es Estabilidad de la cápsula abierta 14 días (es)
- sv Hållbarhet för öppnad flaska 14 dagar (sv)
- tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Sredno
- Sredstvo
- Průměr
- Middelværdi
- Keskmine
- Moyenne
- Mittelwert
- Μέση
- Media
- Nozīmē
- Vidurkis
- Middelverdi
- Srednia
- Média
- Medie
- Среднее
- Srednja vrednost
- Priemer
- Povprečno
- Media
- Medelvärde
- Ortalama

- 3
- Обхват
- Raspon
- Rozsah
- Område
- Vahemik
- Intervalle
- Bereich
- Αναμενόμενο εύρος α
- Intervallo
- Diapazons
- Intervalas
- Område
- Zakres
- Intervalo
- Interval
- диапазон
- Opseg
- Rozsah
- Območje
- Intervalo
- Intervall
- Aralığı



[http://www.medonic.se/
support/downloads](http://www.medonic.se/support/downloads)



Boule

Boule Medical AB
Domnarvsgatan 4
SE-163 53 Spånga, Sweden

1504054

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series

Lot: 22504-61

Prod.: Con-Diff Low

ID: 2250461+

Exp.: 28-Aug-2025

Instr.: M-series

Lot: 22504-62

Prod.: Con-Diff Norm

ID: 2250462+

Exp.: 28-Aug-2025

Instr.: M-series

Lot: 22504-63

Prod.: Con-Diff High

ID: 2250463+

Exp.: 1-Sep-2025

1:



2:



3:



4:



5:



6:



7:



8:



9:

