

Boule Con-5Diff A1

CONTROL

IVD

LOT



Tri-Level	72505-1K	2025/09/05
Low	72505-11	2025/09/05
Normal	72505-12	2025/09/05
High	72505-13	2025/09/05

Open vial stability: 14 days ¹

Medonic M51 / Swelab Lumi

Parameter ²	Unit ³	Low		Normal		High	
		Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵
WBC	10 ⁹ /L	3.71	0.50	8.93	1.00	19.76	2.50
Neu%	%	52.1	16.0	55.6	8.0	61.6	8.0
Lym%	%	35.6	21.0	29.8	8.0	19.5	12.0
Mon%	%	7.0	5.0	7.1	5.0	9.0	6.0
Eos%	%	5.3	5.0	7.5	6.0	9.9	8.0
Bas%	%	5.6	10.0	4.8	10.0	2.5	10.0
Neu#	10 ⁹ /L	1.91	0.50	4.89	0.70	11.98	2.00
Lym#	10 ⁹ /L	1.30	0.70	2.62	0.80	3.78	2.00
Mon#	10 ⁹ /L	0.26	0.30	0.63	0.50	1.75	1.10
Eos#	10 ⁹ /L	0.19	0.20	0.66	0.60	1.92	1.40
Bas#	10 ⁹ /L	0.20	0.40	0.42	0.80	0.48	1.80
RBC	10 ¹² /L	2.33	0.24	4.63	0.30	5.32	0.50
HGB	g/L	62	6	137	6	172	8
HCT	%	19.2	2.0	42.3	3.0	51.5	4.0
MCV	fL	82.5	5.0	91.4	5.0	96.8	6.0
MCH	pg	26.4	2.5	29.6	2.5	32.3	2.5
MCHC	g/L	320	30	324	30	334	30
RDW-CV	%	16.0	3.0	14.7	3.0	14.5	3.0
PLT	10 ⁹ /L	60	20	251	40	491	60
MPV	fL	8.0	3.0	7.9	3.0	8.0	3.0

No mean values and ranges are provided for research parameters-RDW-SD, PDW, PCT, P-LCR, and P-LCC. User may establish these based on their usage.



1	2	3	4	5
bg Трайност на отворен флакон: 14 дни	Параметър	Единица	Средно	Обхват
cs Stabilita otevřené ampulky: 14 dní	Parametr	Jednotka	Průměr	Rozsah
da Åben rør stabilitet: 14 dage	Parameter	Enhed	Middelværdi	Område
de Stabilität geöffneter Flaschen: 14 tage	Parameter	Einheit	Mittelwert	Bereich
el Διατήρηση δείγματος μετά την αποσφράγιση: 14 ημέρα	Παράμετρος	Μονάδα	Μέση	Αναμενόμενο εύρος α
es Estabilidad de la cápsula abierta: 14 días	Parámetro	Unidad	Media	Intervalo
et Stabiilsus avatud vialis: 14 päevade arv	Parameeter	Ühik	Keskmine	Vahemik
fr Stabilité en flacon ouvert: 14 jour	Paramètre	Unité	Moyenne	Intervalle
hr Stabilnost otvorene bočice: 14 dani	Parametar	Jedinica	Sredstvo	Raspon
it Stabilità della fiala aperta: 14 giorni	Parametro	Unità	Media	Intervallo
lt Stabilitumas atidarius buteliuko: 14 dienų	Parametras	Vienetas	Vidurkis	Intervalas
lv Atvērt flakonu stabilitāti: 14 dienas	Parametrs	Vienība	Nozīmē	Diapazons
no Åpen rør/lass stabilitet: 14 dager	Parameter	Enhhet	Middelverdi	Område
pl Trwałość po otwarciu fioleki: 14 Liczba dni	Kryterium	Jednostka	Srednia	Zakres
pt Estabilidade após abertura do frasco: 14 dias	Parâmetro	Unidade	Média	Intervalo
ro Termenul de valabilitate al fiolei destăcut: 14 zile	Parametru	Unitate	Medie	Interval
ru Стабильность вскрытого флакона: 14 дней	Параметр	Единица измерения	Среднее	диапазон
sk Stabilita otvorenej injekčnej liekovky: 14 dní	Parameter	Jednotka	Priemer	Rozsah
sl Stabilnost odprte stekleničke: 14 dni	Parameter	Enota	Povprečno	Območje
sr Proizvod stabilan nakon otvaranja pakovanja: 14 dana	Parametar	Jedinica	Srednja vrednost	Opseg
sv Hållbarhet för öppnad flaska: 14 dagar	Parameter	Enhhet	Medelvärde	Intervall
tr Açık sise stabilitesi gün: 14 gün	Parametre	Birim	Ortalama	Aralığı

1504546 31340-12



www.boule.com

Boule
Boule Medical AB
Domnarvsgatan 4
SE-105 53 Spånga, Sweden