

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	--	--
Low	22506-41	2025-10-31
Normal	22506-42	2025-10-31
High	22506-43	2025-11-04

MEDONIC

CONTROL		L	
LOT		22506-41	
Mean ²		± Range ³	
RBC	10 ¹² /L	2.24	0.12
MCV	fL	66.8	5.0
RDW	%	14.9	4.0
Hct	%	15.0	2.5
Plt	10 ⁹ /L	83	15
MPV	fL	10.7	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	24.6	1.5
MCHC	g/dL	36.7	3.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	44.2	7.0
Mid	%	5.7	5.0
Gran	%	50.1	7.0

CONTROL		N	
LOT		22506-42	
Mean ²		± Range ³	
		4.15	0.18
		79.8	5.0
		14.6	4.0
		33.1	3.0
		238	30
		11.2	1.5
		8.2	0.6
		11.7	0.4
		28.2	1.5
		35.3	3.0
		3.9	0.6
		0.5	0.5
		3.8	1.0
		46.7	5.0
		6.4	5.0
		46.9	5.0

CONTROL		H	
LOT		22506-43	
Mean ²		± Range ³	
		5.06	0.22
		92.7	5.0
		13.8	4.0
		46.9	3.5
		537	60
		11.4	1.5
		20.2	1.8
		16.2	0.5
		32.0	1.5
		34.5	3.0
		12.4	1.5
		1.2	1.0
		6.6	2.5
		61.5	5.0
		5.7	5.0
		32.8	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilitet otvorene bočice 14 dni (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilitumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiolek 14 dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei deschise 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorene injekčneje liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Sredno
 - Sredstvo
 - Průměr
 - Middelverdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervālas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series

Lot: 22506-41

Prod.: Con-Diff Low

ID: 2250641+

Exp.: 31-Oct-2025

Instr.: M-series

Lot: 22506-42

Prod.: Con-Diff Norm

ID: 2250642+

Exp.: 31-Oct-2025

Instr.: M-series

Lot: 22506-43

Prod.: Con-Diff High

ID: 2250643+

Exp.: 4-Nov-2025

1:



B o u l e C B A V 2



B o u l e C B A V 2



B o u l e C B A V 2

2:



2250641+



2250642+



2250643+

3:



25103100001153030000



25103100002153030000



25110400003153030000

4:



10004030100000000000



10004029800000000000



10004029800000000000

5:



22412668500550315025



41518798501170433130



50622927501620546935

6:



08315107152461536730



23830112152821535330



53760114153201534530

7:



03403015040020301708



08206039060050503810



20218124150121006625

8:



14940442700575050170



14640467500645046950



13840615500575032850

9:



11114933524405687255



11114479812885525090



11113970922570793276

Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	--	--
Low	22506-41	2025-10-31
Normal	22506-42	2025-10-31
High	22506-43	2025-11-04

Open vial stability 14 days ¹

SWELAB

CONTROL		L	
LOT		22506-41	
		Mean ²	± Range ³
Swelab Alfa/Alfa Plus	10 ¹² /L	2.25	0.12
RBC	fL	68.6	5.0
MCV	%	14.0	4.0
RDW	%	15.4	2.5
Hct	10 ⁹ /L	87	15
Plt	fL	11.3	1.5
MPV	10 ⁹ /L	3.5	0.3
WBC	g/dL	5.6	0.3
Hgb	pg	24.9	1.5
MCH	g/dL	36.4	3.0
MCHC	10 ⁹ /L	1.7	0.4
Lymph#	10 ⁹ /L	0.1	0.3
Mid#	10 ⁹ /L	1.7	0.8
Gran#	%	47.9	7.0
Lymph	%	3.5	5.0
Mid	%	48.6	7.0
Gran			

CONTROL		N	
LOT		22506-42	
		Mean ²	± Range ³
		4.17	0.18
		81.1	5.0
		13.1	4.0
		33.8	3.0
		244	30
		11.7	1.5
		8.3	0.6
		11.7	0.4
		28.1	1.5
		34.6	3.0
		4.2	0.6
		0.3	0.5
		3.8	1.0
		50.3	5.0
		3.8	5.0
		45.9	5.0

CONTROL		H	
LOT		22506-43	
		Mean ²	± Range ³
		5.05	0.22
		94.1	5.0
		12.5	4.0
		47.5	3.5
		549	60
		11.8	1.5
		20.1	1.8
		16.0	0.5
		31.7	1.5
		33.7	3.0
		12.9	1.5
		0.7	1.0
		6.5	2.5
		64.4	5.0
		3.5	5.0
		32.1	5.0

1	
bg	Трайност на отворен флакон 14 дни (bg)
hr	Stabilnost otvorene bočice 14 dni (hr)
cs	Stabilita otevrené ampulky 14 dny (cs)
da	Åben rør stabilitet 14 dage (da)
et	Stabiilsus avatud viaalis 14 päevade (et)
fr	Stabilité en flacon ouvert 14 jour (fr)
de	Stabilität geöffneter Flaschen 14 tage (de)
el	Διαιτήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it	Stabilità della fiala aperta 14 giorni (it)
lv	Atvērt flakonu stabilitāti 14 dienas (lv)
lt	Stabilumas atidarius buteliuką 14 dienas (lt)
no	Åpen rør/glass stabilitet 14 dager (no)
pl	Trwałość po otwarciu fiolek 14 dni (pl)
pt	Estabilidade após abertura do frasco 14 dias (pt)
ro	Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru	Стабильность вскрытого флакона 14 дней (ru)
sr	Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk	Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl	Stabilita otvorene injekčneje liekovky 14 dni (sl)
es	Estabilidad de la cápsula abierta 14 días (es)
sv	Hållbarhet för öppnad flaska 14 dagar (sv)
tr	Açık şişe stabilitesi 14 gün (tr)

2	
	Средно
	Sredstvo
	Průměr
	Middelverdi
	Keskmine
	Moyenne
	Mittelwert
	Μέση
	Media
	Nožimē
	Vidurkis
	Middelverdi
	Srednia
	Média
	Medie
	Среднее
	Srednja vrednost
	Priemer
	Povprečno
	Media
	Medelvärde
	Ortalama

3	
	Обхват
	Raspon
	Rozsah
	Område
	Vahemik
	Intervalle
	Bereich
	Αναμενόμενο εύρος α
	Intervallo
	Diapazons
	Intervalas
	Område
	Zakres
	Intervalo
	Interval
	диапазон
	Opseg
	Rozsah
	Območje
	Intervalo
	Intervall
	Arahğı

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa

Lot: 22506-41

Prod.: Con-Diff Low

ID: 2250641+

Exp.: 31-Oct-2025

Instr.: Swelab Alfa

Lot: 22506-42

Prod.: Con-Diff Norm

ID: 2250642+

Exp.: 31-Oct-2025

Instr.: Swelab Alfa

Lot: 22506-43

Prod.: Con-Diff High

ID: 2250643+

Exp.: 4-Nov-2025

