

Boule Cal-5Diff G2

CAL

IVD



Calibrator

Lot	
62507-44	2025-10-03

Open vial stability:

5 days¹

Quintus

Parameter ²	RBC 10 ¹² /L	MCV fL	Plt 10 ⁹ /L	MPV fL	WBC 10 ⁹ /L	Hgb g/dL	RDW _{cv} %
Mean ³	4.21	83.7	242	9.5	8.7	13.2	14.8
Range ± ⁴	0.10	2.0	20	1.5	0.3	0.3	1.5

	Target X	Target Y
	85	116

1	2	3	4
bg	Трайност на отворен флакон	14 дни	Обхват
hr	Stabilnost otvorene bočice	14 dni	Raspon
cs	Stabilita otevřené ampulky	14 dny	Rozsah
da	Åben rør stabilitet	14 dage	Område
et	Stabiilsus avatud vialis	14 päevade arv	Vahemik
fr	Stabilité en flacon ouvert	14 jour	Intervalle
de	Stabilität geöffneter Flaschen	14 tage	Bereich
el	Διάρκηση δείγματος μετά την αποσφράγιση	14 ημέρα	Αναμενόμενο εύρος α
it	Stabilità della fiala aperta	14 giorni	Intervallo
lv	Atvērt flakonu stabilitāti	14 dienas	Diapazons
lt	Stabilumas atidarius buteliuką	14 d.	Intervalas
no	Åpen rør/glass stabilitet	14 dager	Område
pl	Trwałość po otwarciu fiolek	14 Liczba dni	Zakres
pt	Estabilidade após abertura do frasco	14 dias	Intervalo
ro	Termenul de valabilitate al fiolei desfăcută	14 zile	Interval
ru	Стабильность вскрытого флакона	14 дней	диапазон
sr	Proizvod stabilan nakon otvaranja pakovanja	14 dana	Opseg
sk	Stabilita otvorenej injekčnej liekovky	14 dní	Rozsah
sl	Stabilnost odprte stekleničke	14 dni	Območje
es	Estabilidad de la cápsula abierta	14 días	Intervalo
sv	Hållbarhet för öppnad flaska	14 dagar	Intervall
tr	Açık şişe stabilitesi	14 gün	Aralığı



www.boule.com

Boule
Boule Medical AB
Domnarvsgatan 4
SE-163 53 Spånga, Sweden

1504393 26577-5

For Quintus with barcode reader

1. In the Main Menu go to “*Calibration QC*” / “*Enter QC reference*”.
2. Highlight the first Barcode field, scan in each calibrator/control barcode of the column in sequence (1-9) and click “OK”.
3. Repeat steps 1 and 2 for each lot to be entered.

Instr.: Quintus

Lot: 62507-44

Prod.: Calibrator

ID: 6250744+

Exp.: 3-Oct-2025

