

LOT: 12505-52
EXP: 2025-10-19
Open Vial Stability: 16 days

EXIGO		Mean	± Range	Expected Range	
WBC	$10^9/L$	8.5	0.6	7.9	9.1
Hct	%	31.0	2.5	28.5	33.5
MCV	fL	74.4	5.0	69.4	79.4
RBC	$10^{12}/L$	4.17	0.18	3.99	4.35
Hgb	g/dL	10.9	0.4	10.5	11.3
MCH	pg	26.1	1.5	24.6	27.6
MCHC	g/dL	35.2	2.2	33.0	37.4
PLT	$10^9/L$	210	30	180	240

This assay sheet is also available on :
www.boule.com

For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12505-52
Prod.: Vet-Con Norm
ID: 1250552+
Exp.: 19-Oct-2025

