

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22508-1K	2026-01-05
Low	22508-11	2026-01-05
Normal	22508-12	2026-01-05
High	22508-13	2026-01-07

MEDONIC

CONTROL		L	
LOT		22508-11	
Medonic M-series/M32		Mean ²	± Range ³
RBC	10 ¹² /L	2.28	0.12
MCV	fL	68.1	5.0
RDW	%	15.4	4.0
Hct	%	15.5	2.5
Plt	10 ⁹ /L	83	15
MPV	fL	11.4	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.8	0.3
MCH	pg	25.4	1.5
MCHC	g/dL	37.4	3.0
Lymph#	10 ⁹ /L	1.4	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.8	0.8
Lymph	%	42.5	7.0
Mid	%	6.2	5.0
Gran	%	51.3	7.0

CONTROL		N	
LOT		22508-12	
		Mean ²	± Range ³
		4.12	0.18
		82.8	5.0
		13.0	4.0
		34.1	3.0
		241	30
		11.6	1.5
		8.3	0.6
		12.0	0.4
		29.1	1.5
		35.2	3.0
		3.8	0.6
		0.6	0.5
		3.9	1.0
		45.9	5.0
		6.9	5.0
		47.2	5.0

CONTROL		H	
LOT		22508-13	
		Mean ²	± Range ³
		5.03	0.22
		91.6	5.0
		13.6	4.0
		46.1	3.5
		539	60
		11.9	1.5
		20.0	1.8
		15.8	0.5
		31.4	1.5
		34.3	3.0
		12.1	1.5
		1.1	1.0
		6.8	2.5
		60.3	5.0
		5.6	5.0
		34.1	5.0

1
bg Трайност на отворен флакон 14 дни (bg)
hr Stabilität otvorene bočice 14 dni (hr)
cs Stabilita otevřené ampulky 14 dní (cs)
da Åben rør stabilitet 14 dage (da)
et Stabilsus avatud viaalis 14 päevade (et)
fr Stabilité en flacon ouvert 14 jour (fr)
de Stabilität geöffneten Flaschen 14 tage (de)
el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it Stabilità della fiala aperta 14 giorni (it)
lv Atvērt flakonu stabilitāti 14 dienas (lv)
lt Stabilitumas atidarius buteliuką 14 dienas (lt)
no Åpen rør/glass stabilitet 14 dager (no)
pl Trwałość po otwarciu fiolek 14 Liczba dni (pl)
pt Estabilidade após abertura do frasco 14 dias (pt)
ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru Стабильность вскрытого флакона 14 дней (ru)
sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl Stabilita otvorene injekčneje liekovky 14 dni (sl)
es Estabilidad de la cápsula abierta 14 días (es)
sv Hållbarhet för öppnad flaska 14 dagar (sv)
tr Açık şişe stabilitesi 14 gün (tr)

2
Средно
Sredstvo
Průměr
Middelværdi
Keskmine
Moyenne
Mittelwert
Μέση
Media
Nozīmē
Vidurkis
Middelverdi
Średnia
Média
Medie
Среднее
Srednja vrednost
Priemer
Povprečno
Media
Medelvärde
Ortalama

3
Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervalo
Intervall
Arahğı

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series

Lot: 22508-11

Prod.: Con-Diff Low

ID: 2250811+

Exp.: 5-Jan-2026

Instr.: M-series

Lot: 22508-12

Prod.: Con-Diff Norm

ID: 2250812+

Exp.: 5-Jan-2026

Instr.: M-series

Lot: 22508-13

Prod.: Con-Diff High

ID: 2250813+

Exp.: 7-Jan-2026

1:



BouleCB AV2



BouleCB AV2



BouleCB AV2

2:



2250811+



2250812+



2250813+

3:



26010500001153030000



26010500002153030000



26010700003153030000

4:



10004029700000000000



10004030100000000000



10004030700000000000

5:



22812681500580315525



41218828501200434130



50322916501580546135

6:



08315114152541537430



24130116152911535230



53960119153141534330

7:



03403014040020301808



08306038060060503910



20018121150111006825

8:



15440425700625051370



13040459500695047250



13640603500565034150

9:



11113296830934995005



11117534496117837296



11112321610884880444

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹

Control	Lot	
Tri-Level	22508-1K	2026-01-05
Low	22508-11	2026-01-05
Normal	22508-12	2026-01-05
High	22508-13	2026-01-07

SWELAB

CONTROL		L	
LOT		22508-11	
Swelab Alfa/Alfa Plus		Mean ²	± Range ³
RBC	10 ¹² /L	2.30	0.12
MCV	fL	69.3	5.0
RDW	%	14.7	4.0
Hct	%	15.9	2.5
Plt	10 ⁹ /L	84	15
MPV	fL	12.1	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.8	0.3
MCH	pg	25.2	1.5
MCHC	g/dL	36.5	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	47.3	7.0
Mid	%	3.7	5.0
Gran	%	49.0	7.0

CONTROL		N	
LOT		22508-12	
		Mean ²	± Range ³
		4.10	0.18
		83.8	5.0
		11.7	4.0
		34.4	3.0
		245	30
		12.3	1.5
		8.2	0.6
		12.0	0.4
		29.3	1.5
		34.9	3.0
		4.1	0.6
		0.3	0.5
		3.8	1.0
		49.8	5.0
		3.9	5.0
		46.3	5.0

CONTROL		H	
LOT		22508-13	
		Mean ²	± Range ³
		5.03	0.22
		93.1	5.0
		12.4	4.0
		46.8	3.5
		557	60
		12.7	1.5
		19.9	1.8
		15.8	0.5
		31.4	1.5
		33.8	3.0
		12.6	1.5
		0.6	1.0
		6.7	2.5
		63.2	5.0
		3.0	5.0
		33.8	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
- hr Stabilitet otvorene bočice 14 dni (hr)
- cs Stabilita otevřených ampulky 14 dní (cs)
- da Åben rør stabilitet 14 dage (da)
- et Stabiilsus avatud vialis 14 päevade (et)
- fr Stabilité en flacon ouvert 14 jour (fr)
- de Stabilität geöffneter Flaschen 14 tage (de)
- el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
- it Stabilità della fiala aperta 14 giorni (it)
- lv Atvērt flakonu stabilitāti 14 dienas (lv)
- lt Stabilitumas atidarius buteliuką 14 dienas (lt)
- no Åpen rør/glass stabilitet 14 dager (no)
- pl Trwałość po otwarciu fiolek 14 dni (pl)
- pt Estabilidade após abertura do frasco 14 dias (pt)
- ro Termenul de valabilitate al fiolei deschise 14 zile (ro)
- ru Стабильность вскрытого флакона 14 дней (ru)
- sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
- sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
- sl Stabilita otvorene injekčneje liekovky 14 dni (sl)
- es Estabilidad de la cápsula abierta 14 días (es)
- sv Hållbarhet för öppnad flaska 14 dagar (sv)
- tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
- Sredstvo
- Průměr
- Middelverdi
- Keskmine
- Moyenne
- Mittelwert
- Μέση
- Media
- Nožimě
- Vidurkis
- Middelverdi
- Średnia
- Média
- Medie
- Среднее
- Srednja vrednost
- Priemer
- Povprečno
- Media
- Medelvärde
- Ortalama

- 3
- Обхват
- Raspon
- Rozsah
- Område
- Vahemik
- Intervalle
- Bereich
- Αναμενόμενο εύρος α
- Intervallo
- Diapazons
- Intervalas
- Område
- Zakres
- Intervalo
- Interval
- диапазон
- Opseg
- Rozsah
- Območje
- Intervalo
- Intervall
- Aralığı

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa

Lot: 22508-11

Prod.: Con-Diff Low

ID: 2250811+

Exp.: 5-Jan-2026

Instr.: Swelab Alfa

Lot: 22508-12

Prod.: Con-Diff Norm

ID: 2250812+

Exp.: 5-Jan-2026

Instr.: Swelab Alfa

Lot: 22508-13

Prod.: Con-Diff High

ID: 2250813+

Exp.: 7-Jan-2026

