

LOT: 12508-52
EXP: 2026-01-18
Open Vial Stability: 16 days

EXIGO		Mean	± Range	Expected Range	
WBC	$10^9/L$	8.6	0.6	8.0	9.2
Hct	%	31.2	2.5	28.7	33.7
MCV	fL	75.2	5.0	70.2	80.2
RBC	$10^{12}/L$	4.15	0.18	3.97	4.33
Hgb	g/dL	10.8	0.4	10.4	11.2
MCH	pg	26.0	1.5	24.5	27.5
MCHC	g/dL	34.6	2.2	32.4	36.8
PLT	$10^9/L$	225	30	195	255

This assay sheet is also available on :
www.boule.com

For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12508-52
Prod.: Vet-Con Norm
ID: 1250852+
Exp.: 18-Jan-2026

