

Boule Con-5Diff A1

CONTROL

IVD

LOT



Tri-Level	72509-1K	2026/01/05
Low	72509-11	2026/01/05
Normal	72509-12	2026/01/05
High	72509-13	2026/01/05

Open vial stability: 14 days ¹

Medonic M51 / Swelab Lumi

Parameter ²	Unit ³	Low		Normal		High	
		Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵
WBC	10 ⁹ /L	3.46	0.50	8.17	1.00	19.14	2.50
Neu%	%	52.0	16.0	57.8	8.0	65.8	8.0
Lym%	%	38.4	21.0	31.0	8.0	20.5	12.0
Mon%	%	6.7	5.0	5.3	5.0	5.8	6.0
Eos%	%	2.9	5.0	5.9	6.0	7.9	8.0
Bas%	%	4.0	10.0	1.0	10.0	2.1	10.0
Neu#	10 ⁹ /L	1.80	0.50	4.72	0.70	12.59	2.00
Lym#	10 ⁹ /L	1.33	0.70	2.53	0.80	3.92	2.00
Mon#	10 ⁹ /L	0.23	0.30	0.43	0.50	1.11	1.10
Eos#	10 ⁹ /L	0.10	0.20	0.48	0.60	1.52	1.40
Bas#	10 ⁹ /L	0.14	0.40	0.08	0.80	0.41	1.80
RBC	10 ¹² /L	2.26	0.24	4.67	0.30	5.39	0.50
HGB	g/L	61	6	136	6	173	8
HCT	%	18.7	2.0	41.5	3.0	51.6	4.0
MCV	fL	82.6	5.0	88.9	5.0	95.8	6.0
MCH	pg	27.0	2.5	29.1	2.5	32.1	2.5
MCHC	g/L	327	30	328	30	335	30
RDW-CV	%	16.0	3.0	15.5	3.0	14.8	3.0
PLT	10 ⁹ /L	63	20	248	40	456	60
MPV	fL	8.3	3.0	8.3	3.0	8.5	3.0

No mean values and ranges are provided for research parameters-RDW-SD, PDW, PCT, P-LCR, and P-LCC. User may establish these based on their usage.



www.boule.com

Support/downloads



1	2	3	4	5
bg	Трайност на отворен флакон: 14 дни	Параметър	Средно	Обхват
cs	Stabilita otevřené ampulky: 14 dny	Parametr	Průměr	Rozsah
da	Åben rør stabilitet: 14 dage	Parameter	Middelværdi	Område
de	Stabilität geöffneter Flaschen: 14 tage	Parameter	Mittelwert	Bereich
el	Διάρθρωση δείγματος μετά την αποσφράγιση: 14 ημέρα	Παράμετρος	Μέση	Αναμενόμενο εύρος α
es	Estabilidad de la cápsula abierta: 14 días	Parámetro	Media	Intervalo
et	Stabiilsus avatud vialis: 14 päevade arv	Parameeter	Keskmine	Vahemik
fr	Stabilité en flacon ouvert: 14 jour	Paramètre	Moyenne	Intervalle
hr	Stabilnost otvorene bočice: 14 dani	Parametar	Sredstvo	Raspon
it	Stabilità della fiala aperta: 14 giorni	Parametro	Media	Intervallo
lt	Stabilumas atidarius buteliuką: 14 dienų	Parametras	Vienetas	Intervalas
lv	Atvērt flakonu stabilitāti: 14 dienas	Parametrs	Vienība	Diapazons
no	Åpen rør/glass stabilitet: 14 dager	Parameter	Enhhet	Område
pl	Trwałość po otwarzeniu fiołki: 14 Liczba dni	Kryterium	Jednostka	Zakres
pt	Estabilidade após abertura do frasco: 14 dias	Parâmetro	Unidade	Intervalo
ro	Termenul de valabilitate al fiolei desfăcute: 14 zile	Parametru	Unitate	Interval
ru	Стабильность вскрытого флакона: 14 дней	Параметр	Единица измерения	диапазон
sk	Stabilita otvorenej injekčnej liekovky: 14 dní	Parameter	Jednotka	Rozsah
sl	Stabilnost odprte stekleničke: 14 dni	Parameter	Enota	Območje
sr	Proizvod stabilan nakon otvaranja pakovanja: 14 dana	Parametar	Jedinica	Opseg
sv	Hållbarhet för öppnad flaska: 14 dagar	Parameter	Enhhet	Intervall
tr	Açık şişe stabilitesi gün: 14 gün	Parametre	Birim	Aralığı