

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹



Control	Lot	
Tri-Level	32509-5K	2026-01-29
Low	32509-51	2026-01-29
Normal	32509-52	2026-01-30
High	32509-53	2026-02-06

MEDONIC

CONTROL		L	
LOT		32509-51	
Medonic M-series/M32		Mean ²	± Range ³
RBC	10 ¹² /L	2.27	0.12
MCV	fL	67.2	7.0
RDW	%	16.2	4.5
Hct	%	15.3	3.0
Plt	10 ⁹ /L	78	15
MPV	fL	11.4	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.7	0.3
MCH	pg	25.1	2.0
MCHC	g/dL	37.3	4.0
Lymph#	10 ⁹ /L	1.4	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	41.8	7.0
Mid	%	6.5	5.0
Gran	%	51.7	7.0

CONTROL		N	
LOT		32509-52	
		Mean ²	± Range ³
		4.19	0.18
		82.2	7.0
		12.9	4.5
		34.4	3.5
		236	30
		11.9	1.5
		8.5	0.6
		12.1	0.4
		28.9	2.0
		35.2	4.0
		3.9	0.6
		0.6	0.5
		4.0	1.0
		45.6	5.0
		7.3	5.0
		47.1	5.0

CONTROL		H	
LOT		32509-53	
		Mean ²	± Range ³
		5.09	0.22
		91.9	7.0
		13.1	4.5
		46.8	4.0
		522	60
		12.1	1.5
		20.4	1.8
		16.1	0.5
		31.6	2.0
		34.4	4.0
		12.3	1.5
		1.2	1.0
		6.9	2.5
		60.1	5.0
		6.1	5.0
		33.8	5.0

1
bg Трайност на отворен флакон 14 дни (bg)
hr Stabilitnost otvorene bočice 14 dani (hr)
cs Stabilita otevrené ampulky 14 dny (cs)
da Åben rør stabilitet 14 dage (da)
et Stabilsus avatud vialis 14 päevade (et)
fr Stabilité en flacon ouvert 14 jour (fr)
de Stabilität geöffneter Flaschen 14 tage (de)
el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it Stabilità della fiala aperta 14 giorni (it)
lv Atvērt flakonu stabilitāti 14 dienas (lv)
lt Stabulumas atidarius buteliuką 14 dienas (lt)
no Åpen rør/glass stabilitet 14 dager (no)
pl Trwałość po otwarciu fiolek 14 dni (pl)
pt Estabilidade após abertura do frasco 14 dias (pt)
ro Termenul de valabilitate al fiolei desfăcută 14 zile (ro)
ru Стабильность вскрытого флакона 14 дней (ru)
sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
es Estabilidad de la cápsula abierta 14 días (es)
sv Hållbarhet för öppnad flaska 14 dagar (sv)
tr Açık şişe stabilitesi 14 gün (tr)

2
Средно
Sredstvo
Průměr
Middelverdi
Keskmine
Moyenne
Mittelwert
Μέση
Media
Nozīmē
Vidurkis
Middelverdi
Średnia
Média
Medie
Среднее
Srednja vrednost
Priemer
Povprečno
Media
Medelvärde
Ortalama

3
Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervalo
Intervall
Aralığı

25122-4



[http://www.medonic.se/
support/download](http://www.medonic.se/support/download)

Boule
Boule Medical AB
Domnarvsgatan 4
SE-163 53 Spånga,
Sweden

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series
Lot: 32509-51
Prod.: Con-Diff Low
ID: 3250951+
Exp.: 29-Jan-2026

Instr.: M-series
Lot: 32509-52
Prod.: Con-Diff Norm
ID: 3250952+
Exp.: 30-Jan-2026

Instr.: M-series
Lot: 32509-53
Prod.: Con-Diff High
ID: 3250953+
Exp.: 6-Feb-2026



3.

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹



Control	Lot	
Tri-Level	32509-5K	2026-01-29
Low	32509-51	2026-01-29
Normal	32509-52	2026-01-30
High	32509-53	2026-02-06

SWELAB

CONTROL		L	
LOT		32509-51	
Swelab Alfa/Alfa Plus		Mean ²	± Range ³
RBC	10 ¹² /L	2.31	0.12
MCV	fL	68.4	7.0
RDW	%	15.2	4.5
Hct	%	15.8	3.0
Plt	10 ⁹ /L	77	15
MPV	fL	12.3	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.7	0.3
MCH	pg	24.7	2.0
MCHC	g/dL	36.1	4.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	46.6	7.0
Mid	%	4.3	5.0
Gran	%	49.1	7.0

CONTROL		N	
LOT		32509-52	
		Mean ²	± Range ³
		4.19	0.18
		83.4	7.0
		11.6	4.5
		34.9	3.5
		243	30
		12.6	1.5
		8.4	0.6
		12.1	0.4
		28.9	2.0
		34.7	4.0
		4.1	0.6
		0.4	0.5
		3.9	1.0
		49.7	5.0
		4.2	5.0
		46.1	5.0

CONTROL		H	
LOT		32509-53	
		Mean ²	± Range ³
		5.07	0.22
		93.0	7.0
		11.8	4.5
		47.2	4.0
		549	60
		12.9	1.5
		20.3	1.8
		16.1	0.5
		31.8	2.0
		34.1	4.0
		12.9	1.5
		0.7	1.0
		6.7	2.5
		63.5	5.0
		3.3	5.0
		33.2	5.0

1
bg Трайност на отворен флакон 14 дни (bg)
hr Stabilitnost otvorene bočice 14 dni (hr)
cs Stabilita otevřené ampulky 14 dní (cs)
da Åben rør stabilitet 14 dage (da)
et Stabiilsus avatud vialis 14 päevade (et)
fr Stabilité en flacon ouvert 14 jour (fr)
de Stabilität geöffneter Flaschen 14 tage (de)
el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it Stabilità della fiala aperta 14 giorni (it)
lv Atvērt flakonu stabilitāti 14 dienas (lv)
lt Stabilitumas atidarius buteliuką 14 dienas (lt)
no Åpen rør/glass stabilitet 14 dager (no)
pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
pt Estabilidade após abertura do frasco 14 dias (pt)
ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru Стабильность вскрытого флакона 14 дней (ru)
sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
es Estabilidad de la cápsula abierta 14 días (es)
sv Hållbarhet för öppnad flaska 14 dagar (sv)
tr Açık şişe stabilitesi 14 gün (tr)

2
Средно
Sredstvo
Průměr
Middelværdi
Keskmine
Moyenne
Mittelwert
Μέση
Media
Nozīmē
Vidurkis
Middelverdi
Średnia
Média
Medie
Среднее
Srednja vrednost
Priemer
Povprečno
Media
Medelvärde
Ortalama

3
Обхват
Raspon
Rozsah
Område
Vahemik
Intervalle
Bereich
Αναμενόμενο εύρος α
Intervallo
Diapazons
Intervalas
Område
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Območje
Intervalo
Intervall
Aralığı

25122-4



[www.swelab.com/
support/downloads](http://www.swelab.com/support/downloads)


Boule Medical AB
Domnarvsgatan 4
SE-163 53 Spånga, Sweden

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa

Lot: 32509-51

Prod.: Con-Diff Low

ID: 3250951+

Exp.: 29-Jan-2026

Instr.: Swelab Alfa

Lot: 32509-52

Prod.: Con-Diff Norm

ID: 3250952+

Exp.: 30-Jan-2026

Instr.: Swelab Alfa

Lot: 32509-53

Prod.: Con-Diff High

ID: 3250953+

Exp.: 6-Feb-2026

1:



B o u l e C B A V 3



B o u l e C B A V 3



B o u l e C B A V 3

2:



3250951+



3250952+



3250953+

3:



26012900001153030000



26013000002153030000



26020600003153030000

4:



20014018000000000000



20014018000000000000



20014018000000000000

5:



23112684700570315830



41918834701210434935



50722930701610547240

6:



07715123152472036140



24330126152892034740



54960129153182034140

7:



03303015040010301708



08406041060040503910



20318129150071006725

8:



15245466700435049170



11645497500425046150



11845635500335033250

9:



11111516114476503617



11119658176906049751



11116636053559085691