

LOT: 12510-15
EXP: 2025-12-24
Open Vial Stability 5 days

EXIGO

	WBC 10 ⁹ /L	RBC 10 ¹² /L	Hgb g/dL	MCV fL	PLT 10 ⁹ /L
System specific values	8.5	4.14	10.4	74.0	227
Tolerance limits ±	0.3	0.10	0.2	2.0	15

This assay sheet is also available on :
www.boule.com

For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo

Lot: 12510-15

Prod.: Vet-Cal

ID: 1251015+

Exp.: 24-Dec-2025

1:



BouleCBAV5

2:



1251015+

3:



25122400034153030000

4:



10004038000000000000

5:



41410740201040299904

6:



22715999069990799908

7:



08503999109991199912

8:



99913999149991599916

9:



11110056527197687659