

# Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	--	--
Low	22512-41	2026-05-04
Normal	22512-42	2026-05-04
High	22512-43	2026-05-04

Open vial stability 14 days <sup>1</sup>

MEDONIC	CONTROL		L
	LOT		22512-41
Medonic M-series/M32		Mean <sup>2</sup>	± Range <sup>3</sup>
RBC	10 <sup>12</sup> /L	2.28	0.12
MCV	fL	69.1	5.0
RDW	%	14.6	4.0
Hct	%	15.8	2.5
Plt	10 <sup>9</sup> /L	79	15
MPV	fL	11.7	1.5
WBC	10 <sup>9</sup> /L	3.2	0.3
Hgb	g/dL	5.8	0.3
MCH	pg	25.4	1.5
MCHC	g/dL	36.7	3.0
Lymph#	10 <sup>9</sup> /L	1.4	0.4
Mid#	10 <sup>9</sup> /L	0.2	0.3
Gran#	10 <sup>9</sup> /L	1.6	0.8
Lymph	%	43.6	7.0
Mid	%	5.7	5.0
Gran	%	50.7	7.0

CONTROL		N
LOT		22512-42
	Mean <sup>2</sup>	± Range <sup>3</sup>
	4.17	0.18
	82.6	5.0
	13.4	4.0
	34.4	3.0
	237	30
	11.2	1.5
	8.1	0.6
	12.0	0.4
	28.8	2.5
	34.9	3.0
	3.7	0.6
	0.5	0.5
	3.9	1.0
	45.5	5.0
	6.2	5.0
	48.3	5.0

CONTROL		H
LOT		22512-43
	Mean <sup>2</sup>	± Range <sup>3</sup>
	5.11	0.22
	92.0	5.0
	13.0	4.0
	47.0	3.5
	535	60
	11.5	1.5
	19.8	1.8
	16.1	0.5
	31.5	1.5
	34.3	3.0
	12.2	1.5
	1.1	1.0
	6.5	2.5
	61.3	5.0
	5.8	5.0
	32.9	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
  - hr Stabilitnost otvorene bočice 14 dani (hr)
  - cs Stabilita otevřené ampulky 14 dní (cs)
  - da Åben rør stabilitet 14 dage (da)
  - et Stabiilsus avatud vialis 14 päevade (et)
  - fr Stabilité en flacon ouvert 14 jour (fr)
  - de Stabilität geöffneter Flaschen 14 tage (de)
  - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
  - it Stabilità della fiala aperta 14 giorni (it)
  - lv Atvērt flakonu stabilitāti 14 dienas (lv)
  - lt Stabilumas atidarius buteliuką 14 dienas (lt)
  - no Åpen rør/glass stabilitet 14 dager (no)
  - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
  - pt Estabilidade após abertura do frasco 14 dias (pt)
  - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
  - ru Стабильность вскрытого флакона 14 дней (ru)
  - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
  - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
  - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
  - es Estabilidad de la cápsula abierta 14 días (es)
  - sv Hållbarhet för öppnad flaska 14 dagar (sv)
  - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
  - Sredstvo
  - Průměr
  - Middelværdi
  - Keskmine
  - Moyenne
  - Mittelwert
  - Μέση
  - Media
  - Nozīmē
  - Vidurkis
  - Middelverdi
  - Średnia
  - Média
  - Medie
  - Среднее
  - Srednja vrednost
  - Priemer
  - Povprečno
  - Media
  - Medelvärde
  - Ortalama

- 3
- Обхват
  - Raspon
  - Rozsah
  - Område
  - Vahemik
  - Intervalle
  - Bereich
  - Αναμενόμενο εύρος α
  - Intervallo
  - Diapazons
  - Intervalas
  - Område
  - Zakres
  - Intervalo
  - Interval
  - диапазон
  - Opseg
  - Rozsah
  - Območje
  - Intervallo
  - Intervall
  - Aralığı

# For Medonic M-series / M32 with barcode reader






**M-Series:**

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

**M32:**

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series	Instr.: M-series	Instr.: M-series
Lot: 22512-41	Lot: 22512-42	Lot: 22512-43
Prod.: Con-Diff Low	Prod.: Con-Diff Norm	Prod.: Con-Diff High
ID: 2251241+	ID: 2251242+	ID: 2251243+
Exp.: 4-May-2026	Exp.: 4-May-2026	Exp.: 4-May-2026

1:  Boul e CBAV2	 Boul e CBAV2	 Boul e CBAV2
2:  2251241+	 2251242+	 2251243+
3:  26050400001153030000	 26050400002153030000	 26050400003153030000
4:  10004030000000000000	 10004029700000000000	 10004030000000000000
5:  22812691500580315825	 41718826501200434430	 51122920501610547035
6:  07915117152541536730	 23730112152881534930	 53560115153151534330
7:  03203014040020301608	 08106037060050503910	 19818122150111006525
8:  14640436700575050770	 13440455500625048350	 13040613500585032950
9:  11118888282012683610	 11116716224863371501	 11110835963494603857

# Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	--	--
Low	22512-41	2026-05-04
Normal	22512-42	2026-05-04
High	22512-43	2026-05-04

Open vial stability 14 days <sup>1</sup>

## SWELAB

Swelab Alfa/Alfa Plus	CONTROL		L
	LOT	22512-41	22512-41
	Mean <sup>2</sup>	± Range <sup>3</sup>	
RBC	10 <sup>12</sup> /L	2.30	0.12
MCV	fL	70.4	5.0
RDW	%	13.9	4.0
Hct	%	16.2	2.5
Plt	10 <sup>9</sup> /L	81	15
MPV	fL	12.5	1.5
WBC	10 <sup>9</sup> /L	3.3	0.3
Hgb	g/dL	5.8	0.3
MCH	pg	25.2	1.5
MCHC	g/dL	35.8	3.0
Lymph#	10 <sup>9</sup> /L	1.6	0.4
Mid#	10 <sup>9</sup> /L	0.1	0.3
Gran#	10 <sup>9</sup> /L	1.6	0.8
Lymph	%	47.7	7.0
Mid	%	4.0	5.0
Gran	%	48.3	7.0

CONTROL	N	
	LOT	22512-42
	Mean <sup>2</sup>	± Range <sup>3</sup>
	4.17	0.18
	83.9	5.0
	12.3	4.0
	35.0	3.0
	243	30
	11.9	1.5
	8.2	0.6
	12.1	0.4
	29.0	1.5
	34.6	3.0
	4.1	0.6
	0.3	0.5
	3.8	1.0
	49.4	5.0
	4.2	5.0
	46.4	5.0

CONTROL	H	
	LOT	22512-43
	Mean <sup>2</sup>	± Range <sup>3</sup>
	5.11	0.22
	93.6	5.0
	12.1	4.0
	47.8	3.5
	550	60
	12.3	1.5
	19.8	1.8
	16.1	0.5
	31.5	1.5
	33.7	3.0
	12.8	1.5
	0.7	1.0
	6.3	2.5
	64.5	5.0
	3.5	5.0
	32.0	5.0

1

bg	Трайност на отворен флакон 14 дни (bg)
hr	Stabilnost otvorene bočice 14 dani (hr)
cs	Stabilita otevřené ampulky 14 dní (cs)
da	Åben rør stabilitet 14 dage (da)
et	Stabiilsus avatud viaalis 14 päevade (et)
fr	Stabilité en flacon ouvert 14 jour (fr)
de	Stabilität geöffneter Flaschen 14 tage (de)
el	Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
it	Stabilità della fiala aperta 14 giorni (it)
lv	Atvērt flakonu stabilitāti 14 dienas (lv)
lt	Stabilumas atidarius buteliuką 14 dienas (lt)
no	Åpen rør/glass stabilitet 14 dager (no)
pl	Trwałość po otwarciu fiołki 14 Liczba dni (pl)
pt	Estabilidade após abertura do frasco 14 dias (pt)
ro	Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
ru	Стабильность вскрытого флакона 14 дней (ru)
sr	Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
sk	Stabilita otvorenej injekčnej liekovky 14 dní (sk)
sl	Stabilita otvorenej injekčnej liekovky 14 dni (sl)
es	Estabilidad de la cápsula abierta 14 días (es)
sv	Hållbarhet för öppnad flaska 14 dagar (sv)
tr	Açık şişe stabilitesi 14 gün (tr)

2

	Средно
	Sredstvo
	Průměr
	Middelverdi
	Keskmine
	Moyenne
	Mittelwert
	Μέση
	Media
	Nozīmē
	Vidurkis
	Middelverdi
	Średnia
	Média
	Medie
	Среднее
	Srednja vrednost
	Priemer
	Povprečno
	Media
	Medelvärde
	Ortalama

3

	Обхват
	Raspon
	Rozsah
	Område
	Vahemik
	Intervalle
	Bereich
	Αναμενόμενο εύρος α
	Intervallo
	Diapazons
	Intervalas
	Område
	Zakres
	Intervalo
	Interval
	диапазон
	Orseg
	Rozsah
	Območje
	Intervalo
	Intervall
	Aralığı

## For Swelab Alfa / Alfa Plus with barcode reader

### Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

### Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa  
 Lot: 22512-41  
 Prod.: Con-Diff Low  
 ID: 2251241+  
 Exp.: 4-May-2026

Instr.: Swelab Alfa  
 Lot: 22512-42  
 Prod.: Con-Diff Norm  
 ID: 2251242+  
 Exp.: 4-May-2026

Instr.: Swelab Alfa  
 Lot: 22512-43  
 Prod.: Con-Diff High  
 ID: 2251243+  
 Exp.: 4-May-2026

