

Boule Con-5Diff A1

CONTROL

IVD

LOT



Tri-Level	72601-1K	2026/05/05
Low	72601-11	2026/05/05
Normal	72601-12	2026/05/05
High	72601-13	2026/05/05

Open vial stability: 14 days ¹

Medonic M51 / Swelab Lumi

Parameter ²	Unit ³	Low		Normal		High	
		Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵	Mean ⁴	± Range ⁵
WBC	10 ⁹ /L	3.44	0.50	8.61	1.00	19.27	2.50
Neu%	%	48.7	16.0	53.8	8.0	61.1	8.0
Lym%	%	39.4	21.0	31.7	8.0	22.2	12.0
Mon%	%	5.7	5.0	6.4	5.0	5.6	6.0
Eos%	%	6.3	5.0	8.2	6.0	11.1	8.0
Bas%	%	0.8	10.0	1.7	10.0	0.8	10.0
Neu#	10 ⁹ /L	1.70	0.50	4.69	0.70	11.91	2.00
Lym#	10 ⁹ /L	1.37	0.70	2.76	0.80	4.34	2.00
Mon#	10 ⁹ /L	0.20	0.30	0.56	0.50	1.09	1.10
Eos#	10 ⁹ /L	0.22	0.20	0.71	0.60	2.17	1.40
Bas#	10 ⁹ /L	0.03	0.40	0.15	0.80	0.15	1.80
RBC	10 ¹² /L	2.33	0.24	4.70	0.30	5.44	0.50
HGB	g/L	61	6	137	6	175	8
HCT	%	19.2	2.0	42.2	3.0	52.6	4.0
MCV	fL	82.4	5.0	89.8	5.0	96.6	6.0
MCH	pg	26.3	2.5	29.1	2.5	32.2	2.5
MCHC	g/L	318	30	325	30	333	30
RDW-CV	%	15.1	3.0	15.0	3.0	14.4	3.0
PLT	10 ⁹ /L	56	20	255	40	477	60
MPV	fL	8.8	3.0	8.6	3.0	8.3	3.0

No mean values and ranges are provided for research parameters-RDW-SD, PDW, PCT, P-LCR, and P-LCC. User may establish these based on their usage.



1
 bg Трайност на отворен флакон: 14 дни
 cs Stabilita otevřené ampulky: 14 dní
 da Åben rør stabilitet: 14 dage
 de Stabilität geöffneter Flaschen: 14 Tage
 el Διατήρηση δείγματος μετά την αποσφράγιση: 14 ημέρες
 es Estabilidad de la cápsula abierta: 14 días
 et Stabiilsus avatud vialis: 14 päevade arv
 fr Stabilité en flacon ouvert: 14 jour
 hr Stabilnost otvorene bočice: 14 dana
 it Stabilità della fiala aperta: 14 giorni
 lt Stabilumas atidarius buteliuką: 14 dienu
 lv Atvērti flakonā stabilitāte: 14 dienas
 no Åpen rør/glass stabilitet: 14 dager
 pl Trwałość po otwarciu fiołki: 14 Liczba dni
 pt Estabilidade após abertura do frasco: 14 dias
 ro Termenul de valabilitate al fiolei desfăcută: 14 zile
 ru Стабильность вскрытого флакона: 14 дней
 sk Stabilita otvorenej injekčnej liekovky: 14 dní
 sl Stabilnost odprte stekleničke: 14 dni
 sr Proizvod stabilan nakon otvaranja pakovanja: 14 dana
 sv Hållbarhet för öppnad flaska: 14 dagar
 tr Açık şişe stabilitesi gün: 14 gün

2
 Параметър
 Parametr
 Parameter
 Parameter
 Παράμετρος
 Parámetro
 Parameeter
 Paramètre
 Parametar
 Parametro
 Parametros
 Parametrs
 Parametr
 Kryterium
 Parámetro
 Parametru
 Параметр
 Parameter
 Parameter
 Parameter
 Parametre

3
 Единица
 Jednotka
 Enhed
 Einheit
 Μονάδα
 Unidad
 Ühik
 Unité
 Jedinica
 Unità
 Vienetas
 Vienība
 Einheit
 Jednostka
 Unidade
 Unitate
 Единица измерения
 Jednotka
 Enota
 Jedinica
 Enhet
 Birim

4
 Средно
 Průměr
 Middelværdi
 Mittelwert
 Μέση
 Media
 Keskmine
 Moyenne
 Sredstvo
 Media
 Vidurkis
 Nozīmē
 Middelverdi
 Średnia
 Média
 Medie
 Среднее
 Priemer
 Povprečno
 Srednja vrednost
 Medelvärde
 Ortalama

5
 Обхват
 Rozsah
 Område
 Bereich
 Αναμενόμενο εύρος α
 Intervalo
 Vahemik
 Intervalle
 Raspon
 Intervallo
 Intervalas
 Diapazons
 Område
 Zakres
 Intervalo
 Intervall
 диапазон
 Rozsah
 Območje
 Opseg
 Intervall
 Aralıqı



www.boule.com

Boule
 Boule Medical AB
 Domnarvsgatan 4
 SE-103 53 Spånga, Sweden

1504546 31340-12

