

LOT: 12602-52
EXP: 2026-07-19
Open Vial Stability: 16 days

EXIGO		Mean	± Range	Expected Range	
WBC	$10^9/L$	8.4	0.6	7.8	9.0
Hct	%	31.7	2.5	29.2	34.2
MCV	fL	75.4	5.0	70.4	80.4
RBC	$10^{12}/L$	4.21	0.18	4.03	4.39
Hgb	g/dL	10.7	0.4	10.3	11.1
MCH	pg	25.4	1.5	23.9	26.9
MCHC	g/dL	33.8	2.2	31.6	36.0
PLT	$10^9/L$	226	30	196	256

This assay sheet is also available on :

www.boule.com

For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo
Lot: 12602-52
Prod.: Vet-Con Norm
ID: 1260252+
Exp.: 19-Jul-2026

