


Boule Cal-5Diff G2

CAL

IVD

Calibrator

Lot	
62604-44	2026-07-03

Open vial stability:

5 days¹

Parameter ²	RBC 10 ¹² /L	MCV fL	Plt 10 ⁹ /L	MPV fL	WBC 10 ⁹ /L	Hgb g/dL	RDWcv %
Mean ³	4.08	84.8	236	9.1	8.7	12.8	14.3
Range ± ⁴	0.10	2.0	20	1.5	0.3	0.3	1.5

Target X	Target Y
90	112

1	2	3	4
bg	14 дни	Параметър	Средно
hr	14 dni	Parameter	Sredstvo
cs	14 dny	Parameter	Průměr
da	14 dage	Parameter	Middelværdi
et	14 päevade arv	Parameeter	Keskmine
fr	14 jour	Paramètre	Moyenne
de	14 tage	Parameter	Mittelwert
el	14 ημέρα	Παράμετρος	Μέση
it	14 giorni	Parametro	Media
lv	14 dienas	Parametrs	Nozīmē
lt	14 d.	Parametras	Vidurkis
no	14 dager	Parameter	Middelverdi
pl	14 Liczba dni	Kryterium	Srednia
pt	14 dias	Parâmetro	Média
ro	14 zile	Parametru	Medie
ru	14 дней	Параметр	Среднее
sr	14 dana	Parameter	Srednja vrednost
sk	14 dní	Parameter	Priemer
sl	14 dni	Parameter	Povprečno
es	14 días	Parámetro	Media
sv	14 dagar	Parameter	Medelvärde
tr	14 gün	Parameter	Ortalama
			Обхват
			Raspon
			Rozsah
			Område
			Vahemik
			Intervalle
			Bereich
			Αναμενόμενο εύρος α
			Intervallo
			Diapazons
			Intervalas
			Område
			Zakres
			Intervalo
			Interval
			диапазон
			Opseg
			Rozsah
			Območje
			Intervalo
			Intervall
			Aralığı



www.boule.com

 **Boule**



Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

26577-6

For Quintus with barcode reader

1. In the Main Menu go to "Calibration QC" / "Enter QC reference".
2. Highlight the first Barcode field, scan in each calibrator/control barcode of the column in sequence (1-9) and click "OK".
3. Repeat steps 1 and 2 for each lot to be entered.

Instr.: Quintus
Lot: 62604-44
Prod.: Calibrator
ID: 6260444+
Exp.: 3-Jul-2026



4.