

Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22604-1K	2026-08-31
Low	22604-11	2026-08-31
Normal	22604-12	2026-08-31
High	22604-13	2026-08-31

Open vial stability 14 days ¹

MEDONIC		CONTROL	L
Medonic M-series/M32		LOT	22604-11
		Mean ²	± Range ³
RBC	10 ¹² /L	2.23	0.12
MCV	fL	65.8	5.0
RDW	%	14.8	4.0
Hct	%	14.7	2.5
Plt	10 ⁹ /L	83	15
MPV	fL	11.1	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.4	0.3
MCH	pg	24.2	1.5
MCHC	g/dL	36.7	3.0
Lymph#	10 ⁹ /L	1.4	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.8	0.8
Lymph	%	43.2	7.0
Mid	%	5.3	5.0
Gran	%	51.5	7.0

CONTROL	N
LOT	22604-12
Mean ²	± Range ³
4.21	0.18
84.2	5.0
13.2	4.0
35.4	3.0
236	30
11.3	1.5
8.5	0.6
12.4	0.4
29.5	1.5
35.0	3.0
4.0	0.6
0.5	0.5
4.0	1.0
46.5	5.0
6.1	5.0
47.4	5.0

CONTROL	H
LOT	22604-13
Mean ²	± Range ³
5.07	0.22
92.4	5.0
13.0	4.0
46.8	3.5
541	60
11.7	1.5
20.7	1.8
16.0	0.5
31.6	1.5
34.2	3.0
12.7	1.5
1.2	1.0
6.8	2.5
61.4	5.0
5.6	5.0
33.0	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelværdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervall
 - Aralığı



www.boule.com

 **Boule**

 Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

25121-8

For Medonic M-series / M32 with barcode reader




























M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series	Instr.: M-series	Instr.: M-series
Lot: 22604-11	Lot: 22604-12	Lot: 22604-13
Prod.: Con-Diff Low	Prod.: Con-Diff Norm	Prod.: Con-Diff High
ID: 2260411+	ID: 2260412+	ID: 2260413+
Exp.: 31-Aug-2026	Exp.: 31-Aug-2026	Exp.: 31-Aug-2026

1:  Boul e CBAV2	 Boul e CBAV2	 Boul e CBAV2
2:  2260411+	 2260412+	 2260413+
3:  26083100001153030000	 26083100002153030000	 26083100003153030000
4:  10004030100000000000	 10004029400000000000	 10004029500000000000
5:  22312658500540314725	 42118842501240435430	 50722924501600546835
6:  08315111152421536730	 23630113152951535030	 54160117153161534230
7:  03403014040020301808	 08506040060050504010	 20718127150121006825
8:  14840432700535051570	 13240465500615047450	 13040614500565033050
9:  11111179484205214713	 11113155227139343615	 11111711821309891331

Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22604-1K	2026-08-31
Low	22604-11	2026-08-31
Normal	22604-12	2026-08-31
High	22604-13	2026-08-31

Open vial stability 14 days ¹

SWELAB

CONTROL		L	
LOT		22604-11	
Swelab Alfa/Alfa Plus		Mean ²	± Range ³
RBC	10 ¹² /L	2.23	0.12
MCV	fL	67.1	5.0
RDW	%	14.0	4.0
Hct	%	15.0	2.5
Plt	10 ⁹ /L	81	15
MPV	fL	11.6	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	24.7	1.5
MCHC	g/dL	36.7	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	46.7	7.0
Mid	%	3.5	5.0
Gran	%	49.8	7.0

CONTROL		N	
LOT		22604-12	
		Mean ²	± Range ³
		4.20	0.18
		85.1	5.0
		12.1	4.0
		35.7	3.0
		238	30
		11.8	1.5
		8.4	0.6
		12.3	0.4
		29.3	1.5
		34.5	3.0
		4.2	0.6
		0.4	0.5
		3.8	1.0
		50.3	5.0
		4.2	5.0
		45.5	5.0

CONTROL		H	
LOT		22604-13	
		Mean ²	± Range ³
		5.00	0.22
		93.9	5.0
		12.0	4.0
		47.0	3.5
		552	60
		12.3	1.5
		20.6	1.8
		15.9	0.5
		31.8	1.5
		33.8	3.0
		13.2	1.5
		0.8	1.0
		6.6	2.5
		64.5	5.0
		3.7	5.0
		31.8	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelværdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Rasjon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Orseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı



www.boule.com

 **Boule**



Boule Medical AB
Fagerstagatan 7
SE-163 53 Spånga,
Sweden

25121-8

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa
Lot: 22604-11
Prod.: Con-Diff Low
ID: 2260411+
Exp.: 31-Aug-2026

Instr.: Swelab Alfa
Lot: 22604-12
Prod.: Con-Diff Norm
ID: 2260412+
Exp.: 31-Aug-2026

Instr.: Swelab Alfa
Lot: 22604-13
Prod.: Con-Diff High
ID: 2260413+
Exp.: 31-Aug-2026

